

Hepsin (M-55): sc-33542

BACKGROUND

Extracellular proteases mediate the digestion of neighboring extracellular matrix components in initial tumor growth, allow desquamation of tumor cells into the surrounding environment, provide the basis for invasion of basement membranes in targeted metastatic organs and are required for release and activation of many growth and angiogenic factors. Hepsin (also known as TMPRSS1) is a type II transmembrane serine protease in mammalian cells that is highly expressed on the surface of hepatocytes. Hepsin is frequently overexpressed in several tumors, suggesting that it is a candidate protease in the invasive process and growth capacity of tumor cells. The basal promoter region of the Hepsin gene contains potential binding sites for SP1, AP2, C/EBP, LF-A1 and E box, which may be responsible for the ubiquitous expression of the protein, which shows preferential expression in liver and kidney. Hepsin is located at the plasma membrane, with its catalytic subunit (C-terminal half) at the cell surface and its N-terminus facing the cytosol. Hepsin has been shown to play a role in normal cell growth, embryogenesis, hepatocyte growth, blood coagulation and fertilization. In addition, Hepsin converts zymogen Factor VII to Factor VIIa, which is capable of initiating a coagulation pathway on the cell surface and ultimately leads to Thrombin formation.

REFERENCES

1. Tsuji, A., et al. 1991. Characterization of Hepsin, a membrane bound protease. *Biomed. Biochim. Acta* 50: 791-793.
2. Tsuji, A., et al. 1991. Hepsin, a cell membrane-associated protease. Characterization, tissue distribution, and gene localization. *J. Biol. Chem.* 266: 16948-16953.
3. Kazama, Y., et al. 1995. Hepsin, a putative membrane-associated serine protease, activates human factor VII and initiates a pathway of blood coagulation on the cell surface leading to thrombin formation. *J. Biol. Chem.* 270: 66-72.
4. Tanimoto, H., et al. 1997. Hepsin, a cell surface serine protease identified in hepatoma cells, is overexpressed in ovarian cancer. *Cancer Res.* 57: 2884-2887.

CHROMOSOMAL LOCATION

Genetic locus: HPN (human) mapping to 19q13.12; Hpn (mouse) mapping to 7 B1.

SOURCE

Hepsin (M-55) is a rabbit polyclonal antibody raised against amino acids 221-275 mapping within an internal region of Hepsin of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

Hepsin (M-55) is recommended for detection of Hepsin catalytic chain of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Hepsin (M-55) is also recommended for detection of Hepsin catalytic chain in additional species, including canine, bovine and porcine.

Suitable for use as control antibody for Hepsin siRNA (h): sc-60782, Hepsin siRNA (m): sc-41657, Hepsin shRNA Plasmid (h): sc-60782-SH, Hepsin shRNA Plasmid (m): sc-41657-SH, Hepsin shRNA (h) Lentiviral Particles: sc-60782-V and Hepsin shRNA (m) Lentiviral Particles: sc-41657-V.


Molecular Weight of Hepsin: 51 kDa.

Positive Controls: WEHI-231 whole cell lysate: sc-2213, NIH/3T3 whole cell lysate: sc-2210 or MDCK cell lysate: sc-2252.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



Hepsin (M-55): sc-33542. Western blot analysis of Hepsin expression in NIH/3T3 whole cell lysate.



Hepsin (M-55): sc-33542. Western blot analysis of Hepsin expression in WEHI-231 whole cell lysate.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **Hepsin (2D5): sc-517056**, our highly recommended monoclonal alternative to Hepsin (M-55).