

# HMW-Guanylin (K-14): sc-33833

## BACKGROUND

The family of Guanylin regulatory peptides, including Guanylin and Uroguanylin, are strongly expressed in intestinal mucosa and regulate intestinal fluid secretion during digestion. Guanylins are also involved in acid neutralization and the regulation of membrane-bound Guanylate cyclase signaling molecules. Guanylin and Uroguanylin are secreted primarily in the stomach, intestine and colon. Guanylin is also detected in plasma. Guanylin is an endogenous activator of intestinal guanylate cyclase. It stimulates intestinal guanylate cyclase through the same receptor binding region as the heat-stable enterotoxins. Gut enterochromaffin cells synthesize Guanylin to be a prohormone of 115 amino acids which is then processed to the molecular form of 94 amino acids. This form is found circulating in the blood.

## REFERENCES

- Currie, M.G., Fok, K.F., Kato, J., Moore, R.J., Hamra, F.K., Duffin, K.L. and Smith, C.E. 1992. Guanylin: an endogenous activator of intestinal guanylate cyclase. *Proc. Natl. Acad. Sci. USA* 89: 947-951.
- de Sauvage, F.J., Keshav, S., Kuang, W.J., Gillett, N., Henzel, W. and Goeddel, D.V. 1992. Precursor structure, expression and tissue distribution of human Guanylin. *Proc. Natl. Acad. Sci. USA* 89: 9089-9093.
- Wiegand, R.C., Kato, J. and Currie, M.G. 1992. Rat guanylin cDNA: characterization of the precursor of an endogenous activator of intestinal guanylate cyclase. *Biochem. Biophys. Res. Commun.* 185: 812-817.
- Schulz, S., Chrisman, T.D. and Garbers, D.L. 1992. Cloning and expression of Guanylin. Its existence in various mammalian tissues. *J. Biol. Chem.* 267: 16019-16021.
- Wiegand, R.C., Kato, J., Huang, M.D., Fok, K.F., Kachur, J.F. and Currie, M.G. 1992. Human Guanylin: cDNA isolation, structure and activity. *FEBS Lett.* 311: 150-154.
- Kuhn, M., Raida, M., Adermann, K., Schulz-Knappe, P., Gerzer, R., Heim, J.M. and Forssmann, W.G. 1993. The circulating bioactive form of human Guanylin is a high molecular weight peptide (10.3 kDa). *FEBS J.* 318: 205-209.
- Skelton, N.J., Garcia, K.C., Goeddel, D.V., Quan, C. and Burnier, J.P. 1994. Determination of the solution structure of the peptide hormone Guanylin: observation of a novel form of topological stereoisomerism. *Biochemistry* 33: 13581-13592.
- Hill, O., Kuhn, M., Zucht, H.D., Cetin, Y., Kulaksiz, H., Adermann, K., Klock, G., Rechkemmer, G., Forssmann, W.G. and Magert, H.J. 1995. Analysis of the human Guanylin gene and the processing and cellular localization of the peptide. *Proc. Natl. Acad. Sci. USA* 92: 2046-2050.

## CHROMOSOMAL LOCATION

Genetic locus: GUCA2A (human) mapping to 1p35-p34

## SOURCE

HMW-Guanylin (K-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of HMW-Guanylin of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-33833 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

HMW-Guanylin (K-14) is recommended for detection of HMW-Guanylin of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for HMW-Guanylin siRNA (h): sc-45332, HMW-Guanylin shRNA Plasmid (h): sc-45332-SH and HMW-Guanylin shRNA (h) Lentiviral Particles: sc-45332-V.

Molecular Weight of HMW-Guanylin: 13 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotting A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.


 MONOS  
Satisfaction  
Guaranteed

Try **HMW-Guanylin (MAK-L-G-11): sc-59559**, our highly recommended monoclonal alternative to HMW-Guanylin (K-14).