

## VP16 (H-15): sc-34070

### BACKGROUND

The GAL4 protein of *Saccharomyces cerevisiae* is one of the most thoroughly characterized transcriptional activators. Since the N-terminal 147 amino acid residues of GAL4 are sufficient to mediate specific and strong binding to DNA, but are incapable of efficient transcriptional activation, this protein fragment has frequently been used to confer specific DNA binding in experiments examining transcriptional activation functions of heterologous proteins. This approach is facilitated by the finding that higher eukaryotes lack endogenous proteins that enhance transcription from the consensus GAL4-binding site. Fusions between GAL4 (amino acids 1-147) and activating domains from a variety of transcriptional regulatory proteins can activate transcription in yeast, plant, insects and mammalian cells. A unique "two-hybrid" system has been developed using GAL4 fusions in yeast to identify specific protein-protein interactions. Another "two-hybrid" system utilizes the DNA binding domain of the *E. coli* protein Lex A and the transactivity domain of the HSV protein VP16.

### REFERENCES

1. Johnston, M. 1987. A model fungal gene regulatory mechanism: the GAL genes of *Saccharomyces cerevisiae*. *Microbiol. Rev.* 51: 458-476.
2. Ma, J., et al. 1987. Deletion analysis of GAL4 defines two transcriptional activating segments. *Cell* 48: 847-853.
3. Fields, S., et al. 1989. A novel genetic system to detect protein-protein interactions. *Nature* 340: 245-246.
4. Ptashne, M., et al. 1990. Activators and targets. *Nature* 346: 329-331.
5. Song, O., et al. 1991. Pheromone-dependent phosphorylation of the yeast STE12 protein correlates with transcriptional activation. *Genes Dev.* 5: 741-750.
6. Chien, C., et al. 1991. The two-hybrid system; a method to identify and clone genes for proteins that interact with a protein of interest. *Proc. Natl. Acad. Sci. USA* 88: 9578-9582.
7. Sadowski, I., et al. 1992. GAL4 fusion vectors for expression in yeast or mammalian cells. *Gene* 118: 137-141.
8. Fields, S., et al. 1994. The two-hybrid system: an assay for protein-protein interactions. *Trends Gen.* 10: 286-292.

### SOURCE

VP16 (H-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of VP16 of Herpes Virus origin.

### PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-34070 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### APPLICATIONS

VP16 (H-15) is recommended for detection of VP16 and VP16 fusion proteins of N/A origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.



Try **VP16 (1-21): sc-7545** or **VP16 (14-5): sc-7546**, our highly recommended monoclonal alternatives to VP16 (H-15). Also, for AC, HRP, FITC, PE, Alexa Fluor® 488 and Alexa Fluor® 647 conjugates, see **VP16 (1-21): sc-7545**.