

LXR β (S-16): sc-34343

BACKGROUND

Retinoids are metabolites of vitamin A (retinol) and are believed to represent important signaling molecules during vertebrate development and tissue differentiation. The cooperation of liver X receptors (LXRs) α and β and retinoic X receptor (RXR) modulate the expression of several genes involved in lipid metabolism in hepatocytes and macrophages. RXR is the receptor for 9-*cis* retinoic acid and dimerizes with VDR, TR, PPAR and several novel receptors including liver X receptors LXR α (also referred to as RLD-1), LXR β and FXR. FXR and LXR fall into a category of proteins termed "orphan receptors" because of their lack of a defined function and, in the case of LXR, the lack of a defined ligand. Both LXR/RXR and FXR/RXR heterodimers retain their responsiveness to 9-*cis* retinoic acid. LXR α and LXR β share considerable sequence homology and several functions, respond to the same endogenous and synthetic ligands, and play critical roles in maintaining lipid homeostasis. LXR β is ubiquitously expressed and enriched in tissues of neuronal and endocrine origin.

REFERENCES

1. Mangelsdorf, D.J., et al. 1994. The retinoid receptors. In Sporn, M.B., et al, eds. The Retinoids: Biology, Chemistry, and Medicine. New York: Raven Press, 319-349.
2. Bhat, M.K., et al. 1994. Phosphorylation enhances the target gene sequence-dependent dimerization of thyroid hormone receptor with retinoid X receptor. Proc. Natl. Acad. Sci. USA 91: 7927-7931.

CHROMOSOMAL LOCATION

Genetic locus: NR1H2 (human) mapping to 19q13.33; Nr1h2 (mouse) mapping to 7 B4.

SOURCE

LXR β (S-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of LXR β of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin. Also available as TransCruz reagent for Gel Supershift and ChIP applications, sc-34343 X, 200 μ g/0.1 ml.

Blocking peptide available for competition studies, sc-34343 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

APPLICATIONS

LXR β (S-16) is recommended for detection of LXR β of human and, to a lesser extent, mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


Suitable for use as control antibody for LXR β siRNA (h): sc-45316, LXR β siRNA (m): sc-45317, LXR β shRNA Plasmid (h): sc-45316-SH, LXR β shRNA Plasmid (m): sc-45317-SH, LXR β shRNA (h) Lentiviral Particles: sc-45316-V and LXR β shRNA (m) Lentiviral Particles: sc-45317-V.

LXR β (S-16) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.


Molecular Weight of LXR β : 56 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227, LXR β (h): 293T Lysate: sc-112157 or HeLa nuclear extract: sc-2120.

DATA



LXR β (S-16): sc-34343. Western blot analysis of LXR β expression in non-transfected: sc-117752 (A) and human LXR β transfected: sc-112157 (B) 293T whole cell lysates.



LXR β (S-16): sc-34343. Western blot analysis of LXR β expression in HeLa nuclear extract.

SELECT PRODUCT CITATIONS

1. Shen, Q., et al. 2011. Liver X receptor-retinoid X receptor (LXR-RXR) heterodimer cistrome reveals coordination of LXR and AP1 signaling in keratinocytes. J. Biol. Chem. 286: 14554-14563.



Try **LXR α / β (H-7): sc-377260** or **LXR β (H-8): sc-133221**, our highly recommended monoclonal alternatives to LXR β (S-16). Also, for AC, HRP, FITC, PE, Alexa Fluor[®] 488 and Alexa Fluor[®] 647 conjugates, see **LXR α / β (H-7): sc-377260**.