

# CCL14 (N-13): sc-34876

## BACKGROUND

Chemokines are members of a superfamily of inducible, secreted, pro-inflammatory cytokines. Members of the chemokine family exhibit 20 to 50% homology in their predicted amino acid sequences and are divided into four subfamilies. CCL14 belongs to the intercrine  $\beta$  (chemokine C-C) family. CCL14 has weak activities on human monocytes and acts via receptors that also recognize MIP-1 $\alpha$ . CCL14 also enhances the proliferation of CD34 myeloid progenitor cells. The processed form of CCL14, designated HCC-1(9-74), is a chemotactic factor that attracts monocytes, eosinophils and T cells and is a ligand for CCR1, CCR3 and CCR5. Various membrane-associated and soluble proteases selectively cleave specific chemokines. Precursor plasma chemokines (CXCL7, CCL14) need to be proteolytically processed to obtain receptor affinity.

## REFERENCES

1. Munch, J., et al. 2002. Hemofiltrate CC chemokine 1[9-74] causes effective internalization of CCR5 and is a potent inhibitor of R5-tropic human immunodeficiency virus type 1 strains in primary T cells and macrophages. *Antimicrob. Agents Chemother.* 46: 982-990.
2. Forssmann, U., et al. 2004. n-Nonanoyl-CC chemokine ligand 14, a potent CC chemokine ligand 14 analogue that prevents the recruitment of eosinophils in allergic airway inflammation. *J. Immunol.* 173: 3456-3466.
3. Shen, Y., et al. 2004. Distinct gene expression profiles in different B cell compartments in human peripheral lymphoid organs. *BMC Immunol.* 5: 20.
4. Van Damme, J., et al. 2004. Chemokine-protease interactions in cancer. *Semin. Cancer Biol.* 14: 201-208.
5. Blain, K.Y., et al. 2007. Structural and functional characterization of CC chemokine CCL14. *Biochemistry* 46:10008-10015.

## CHROMOSOMAL LOCATION

Genetic locus: CCL14 (human) mapping to 17q12.

## SOURCE

CCL14 (N-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of CCL14 of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-34876 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.

## APPLICATIONS

CCL14 (N-13) is recommended for detection of CCL14 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for CCL14 siRNA (h): sc-45580, CCL14 shRNA Plasmid (h): sc-45580-SH and CCL14 shRNA (h) Lentiviral Particles: sc-45580-V.

Molecular Weight of CCL14: 9 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.