


C T International

Your Hand and Body Protection Specialist

Product Specifications



LEX Latex Lightly Powdered Exam

Physical Properties

Smooth Surface

Tapered, Beaded Cuff

10" Length, Ambidextrous

Leak Tested

Meets or Exceeds FDA Requirements with 510(k) Registration

Meets or Exceeds ASTM D3578-95 Standards

Sampled in Accordance with MIL STD-105E

Complies with U.S. Government 21CFR 177.2600 for Food Handling

Catalog Numbers and Packaging

Length: 10" (± 0.5)

Thickness: 5.0 mil (± 1)

Tensile Strength: 3400 psi/23 mpa

Elongation: 745%

Packaging: Dispenser Boxed, 100 Gloves/Box, 1000 Gloves/Case

Cat No. /Sizes: LEX100 - X-Small, LEX101 - Small, LEX102 - Medium,
LEX103 - Large, LEX104 - X-Large

Note: All specifications subject to change without notice



C T International

800-755-7575 fax 805-544-5572

4340 Old Santa Fe Road San Luis Obispo, CA 93401