

Products designed with a purpose...
...to help you help them.



CAPIJECT® Capillary Micro Collection Tubes

- Extra-wide, built-in collection lip for enhanced specimen flow
- Single piece, tapered design for easy collection and handling
- Color-coded fill lines and caps for volume accuracy and easy identification

With a full compliment of collection tubes for serum, plasma and whole blood determinations, a small but adequate amount of blood can be easily obtained by skin puncture.

LATEX FREE

 **TERUMO®**

TERUMO MEDICAL PRODUCTS

CAPIJECT® Capillary Micro Collection Tubes

CAPIJECT Makes it Easy

- Extra-wide, built-in collection lip for enhanced specimen flow
- Single piece, tapered design for easy collection and handling
- Color-coded fill lines and caps for volume accuracy and easy identification
- Lot number and expiration date printed on each tube for easy identification
- Silicone coated anticoagulant tubes designed to reduce clotting
- Olefin-based gel for clean separations
- Tube nests in cap allowing tube to stand alone
- Certified for lead testing
- Latex-free



Catalog Number	Tube Color	Tube additive	Application	Lead Testing Certified	Fill volume μL
T-MG	Cherry Red	Gel Barrier with clot activator (silica particles and inert gel)	Serum determinations in clinical chemistry	X	500
T-MQK	Lavender	K2-EDTA (0.78 mg)	Whole blood hematology procedures	X	500
T-M	Brick Red	No Additive/Clot Tube	Serum determinations in clinical chemistry, serology and blood banking	X	500
T-MLH	Green	Lithium Heparin (12.5 IU)	Plasma determinations in clinical chemistry	X	500
T-MLHG	Mint Green	Lithium Heparin with Gel (12.5 IU)	Plasma determinations in clinical chemistry	X	500
T-MLHGA	Mint Green	Amber tube with Lithium Heparin and gel barrier (12.5 IU)	Light-sensitive plasma determinations in clinical chemistry	X	500
T-MPS	Gray	Sodium Fluoride (0.72 mg) Potassium Oxalate (0.58 mg)	Glucose determinations	X	250
T-MGA	Gold	Amber tube with gel barrier and clot activator (silica particles and inert gel)	Light-sensitive serum analyte determinations (Bilirubin determinations)		500

FOR IN VITRO DIAGNOSTIC USE

Packaging

- 50 Capiject Capillary Blood Collection Tubes per bag
- 5 bags per box - 250 Capiject Tubes
- 4 boxes per case - 1,000 Capiject Tubes

TerumoTMP.com/Capiject

US Office

2101 Cottontail Lane
Somerset, NJ 08873
Phone: 732.302.4900
Fax: 732.302.3093
Toll free: 800.283.7866
Fax for Orders: 800.411.5870

Canada Office

Terumo Medical Corporation
C/O UPS - Supply Chain Solutions
1165 Creditstone Road
Vaughan, Ontario L4K 4N7
Phone: 800.668.5489
Fax: 905.760.1215



Ⓢ and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. CAPIJECT is a registered trademark of Terumo Medical Corporation, Somerset, NJ. TERUMO is registered trademark of Terumo Corporation, Tokyo, Japan ©2005 Terumo Medical Corporation