

World Champion Pepi Coat Conditioner

Style Stable Products
Instant Sheen For That Show/Ring Look

An easy to use conditioner and glosser to give that world champion style show ring look. Extra gentle and will not irritate animal's skin.

Directions

Spray coat evenly with sweeping strokes, holding can 12 to 15 inches away. Wipe it with a towel or a brush and re-apply if necessary.

DO NOT spray near eyes or nostrils.

WARNING

FLAMMABLE. DO NOT USE NEAR FIRE OR FLAME OR WHILE SMOKING.
CONTENTS UNDER PRESSURE

Do not puncture or incinerate container. Do not store at temperatures above 120° F (48.89° C). Use only as directed. Intentional misuse by deliberately concentrating and inhaling the contents can be harmful or fatal.

Keep this and all aerosols out of reach of children

Avoid excessive inhalation

Use sparingly in well-ventilated area.

Style Stable Products
11 N Main St
Hutchinson, MN 55350
1-320-587-3217
www.stylestableblanket.com

Disclaimer

Every effort has been made by Jeffers to ensure the accuracy of the information listed above. However, it remains your responsibility to become familiar with the products you are purchasing. Please consult your family veterinarian if you have any questions.