## SANTA CRUZ BIOTECHNOLOGY, INC.

# HEB (A-6): sc-365980



## BACKGROUND

Differentiation of myogenic cells is regulated by multiple positively and negatively acting factors. One well characterized family of helix-loop-helix (HLH) proteins known to play an important role in the regulation of muscle cell development includes Myo D, myogenin, Myf-5 and Myf-6 (also designated MRF-4 or herculin). Myo D transcription factors form heterodimers with products of a more widely expressed family of bHLH genes, the E family, which consists of at least three distinct genes: E2A, IF2 and HEB. Myo D-E heterodimers bind avidly to consensus (CANNTG) E box target sites that are functionally important elements in the upstream regulatory sequences of many muscle-specific terminal differentiation genes. Both homo- and hetero-oligomers of these proteins are able to distinguish very closely related E box proteins and are believed to play important roles in lineage specific gene expression.

# REFERENCES

- 1. Braun, T., et al. 1989. A novel human muscle factor related to but distinct from MyoD1 induces myogenic conversion in 10T1/2 fibroblasts. EMBO J. 8: 701-709.
- 2. Wright, W.E., et al. 1989. Myogenin, a factor regulating myogenesis, has a domain homologous to MyoD. Cell 56: 607-617.
- 3. Murre, C., et al. 1989. Interactions between heterologous helix-loop-helix proteins generate complexes that bind specifically to a common DNA sequence. Cell 58: 537-544.
- 4. Rhodes, S.J., et al. 1989. Identification of MRF4: a new member of the muscle regulatory factor gene family. Genes Dev. 3: 2050-2061.
- 5. Miner, J.H., et al. 1990. Herculin, a fourth member of the MyoD family of myogenic regulatory genes. Proc. Natl. Acad. Sci. USA 87: 1089-1093.

#### **CHROMOSOMAL LOCATION**

Genetic locus: TCF12 (human) mapping to 15g21.3; Tcf12 (mouse) mapping to 9 D.

#### SOURCE

HEB (A-6) is a mouse monoclonal antibody specific for an epitope mapping between amino acids 652-683 within the C-terminus of HEB of human origin.

#### PRODUCT

Each vial contains 200  $\mu$ g IgG<sub>3</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin. Also available as TransCruz reagent for Gel Supershift and ChIP applications, sc-365980 X, 200 µg/0.1 ml.

Blocking peptide available for competition studies, sc-365980 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% stabilizer protein).

#### **STORAGE**

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

#### **APPLICATIONS**

HEB (A-6) is recommended for detection of HEB of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

HEB (A-6) is also recommended for detection of HEB in additional species, including equine and bovine.

Suitable for use as control antibody for HEB siRNA (h): sc-35552, HEB siRNA (m): sc-35553, HEB shRNA Plasmid (h): sc-35552-SH, HEB shRNA Plasmid (m): sc-35553-SH, HEB shRNA (h) Lentiviral Particles: sc-35552-V and HEB shRNA (m) Lentiviral Particles: sc-35553-V.

HEB (A-6) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of HEB: 85 kDa.

Positive Controls: Jurkat nuclear extract: sc-2132, HeLa nuclear extract: sc-2120 or HeLa whole cell lysate: sc-2200.

### **RECOMMENDED SUPPORT REAGENTS**

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG K BP-HRP: sc-516102 or m-IgG K BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

#### DATA





HEB (A-6): sc-365980. Western blot analysis of HEB expression in HeLa (A), Jurkat (B) and Hep G2 (C) nuclear extracts

HEB (A-6): sc-365980. Western blot analysis of HEB expression in HeLa whole cell lysate

#### **SELECT PRODUCT CITATIONS**

1. Wang, L., et al. 2023. ASCL1 characterizes adrenergic neuroblastoma via its pioneer function and cooperation with core regulatory circuit factors. Cell Rep. 42: 113541.

### **RESEARCH USE**

For research use only, not for use in diagnostic procedures.