

Calpain 6 (H-90): sc-366186

BACKGROUND

Calpains are calcium-activated thiol proteases involved in intracellular processing of proteins and signal transduction. The classic calpains are heterodimers with one large subunit, one small subunit and five EF-hand-calcium binding structures. The large subunit varies between family members and can be active without the small subunit. Calpain 6 (Capn6), also designated calpamodulin, CalpM or calpain-like protease X-linked, belongs to the peptidase C-2 family of proteins. It is exclusively expressed in placenta and may play a role in embryogenesis. Calpain 6, in addition to Calpain 10 and Calpain 5, differs from other subfamily members by the presence of a T-domain that is homologous to the *Caenorhabditis elegans* protein Tra-3. Calpain 6 is distinct from other calpains in that it lacks the cysteine and histidine residues required for protease activity. Overexpression of Calpain 6 leads to failure of cytokinesis and formation of microtubule bundles. This results in cells containing multiple nuclei. Calpain 6 may function in Actin organization and microtubule stability.

REFERENCES

- Zhao, P., Iezzi, S., Carver, E., Dressman, D., Gridley, T., Sartorelli, V. and Hoffman, E.P. 2002. SLUG is a novel downstream target of MyoD. Temporal profiling in muscle regeneration. *J. Biol. Chem.* 277: 30091-30101.
- Czogalla, A. and Sikorski, A.F. 2005. Spectrin and Calpain: a "target" and a "sniper" in the pathology of neuronal cells. *Cell. Mol. Life Sci.* 62: 1913-1924.
- Markmann, A., Schäfer, S., Linz, W., Löhn, M., Busch, A.E. and Wohlfart, P. 2005. Downregulation of Calpain 9 is linked to hypertensive heart and kidney disease. *Cell. Physiol. Biochem.* 15: 109-116.
- Leung-Toung, R., Zhao, Y., Li, W., Tam, T.F., Karimian, K. and Spino, M. 2006. Thiol proteases: inhibitors and potential therapeutic targets. *Curr. Med. Chem.* 13: 547-581.
- Paquet-Durand, F., Azadi, S., Hauck, S.M., Ueffing, M., van Veen, T. and Ekström, P. 2006. Calpain is activated in degenerating photoreceptors in the rd1 mouse. *J. Neurochem.* 96: 802-814.
- Ersfeld, K. and Croall, D.E. 2007. The calpains: modular designs and functional diversity. *Genome Biol.* 8: 218-218.
- Tonami, K., Kurihara, Y., Aburatani, H., Uchijima, Y., Asano, T. and Kurihara, H. 2007. Calpain 6 is involved in microtubule stabilization and cytoskeletal organization. *Mol. Cell. Biol.* 27: 2548-2561.
- Lee, S.J., Choi, Y.L., Lee, E.J., Kim, B.G., Bae, D.S., Ahn, G.H. and Lee, J.H. 2007. Increased expression of Calpain 6 in uterine sarcomas and carcinosarcomas: an immunohistochemical analysis. *Int. J. Gynecol. Cancer* 17: 248-253.
- Evans, J.S. and Turner, M.D. 2007. Emerging functions of the calpain superfamily of cysteine proteases in neuroendocrine secretory pathways. *J. Neurochem.* 103: 849-859.

CHROMOSOMAL LOCATION

Genetic locus: CAPN6 (human) mapping to Xq23; Capn6 (mouse) mapping to X F2.

SOURCE

Calpain 6 (H-90) is a rabbit polyclonal antibody raised against amino acids 391-480 mapping within an internal region of Calpain 6 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

Calpain 6 (H-90) is recommended for detection of Calpain 6 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Calpain 6 (H-90) is also recommended for detection of Calpain 6 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for Calpain 6 siRNA (h): sc-62066, Calpain 6 siRNA (m): sc-62067, Calpain 6 shRNA Plasmid (h): sc-62066-SH, Calpain 6 shRNA Plasmid (m): sc-62067-SH, Calpain 6 shRNA (h) Lentiviral Particles: sc-62066-V and Calpain 6 shRNA (m) Lentiviral Particles: sc-62067-V.

Molecular Weight of Calpain 6: 74 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.