

# NKHC2 (E-5): sc-374468

## BACKGROUND

Neuronal kinesin heavy chain 2 (NKHC2) is a 1,032 amino acid protein that is part of the kinesin superfamily which consists of the heavy chains of conventional kinesin. NKHC is expressed throughout the central nervous system, but is highly expressed in certain subsets of neurons. NKHC has a unique C-terminal stretch of 69 amino acids and interacts with dystrobrevin, an adaptor/scaffolding protein. This interaction may play a role in the transport and targeting of components of the dystrophin-associated protein complex to precise sites in the cell. NKHC may also be involved in the microtubule-dependent slow axonal transport of neurofilament proteins during the maturation of neuronal cells.

## REFERENCES

1. Vignali, G., et al. 1997. Expression of neuronal kinesin heavy chain is developmentally regulated in the central nervous system of the rat. *J. Neurochem.* 69: 1840-1849.
2. Rahman, A., et al. 1999. Defective kinesin heavy chain behavior in mouse kinesin light chain mutants. *J. Cell Biol.* 146: 1277-1288.
3. Li, J.Y., et al. 1999. Axonal transport and distribution of immunologically distinct kinesin heavy chains in rat neurons. *J. Neurosci. Res.* 58: 226-241.
4. Kanai, Y., et al. 2000. KIF5C, a novel neuronal kinesin enriched in motor neurons. *J. Neurosci.* 20: 6374-6384.
5. Cai, Y., et al. 2001. The docking of kinesins, KIF5B and KIF5C, to Ran-binding protein 2 (RanBP2) is mediated via a novel RanBP2 domain. *J. Biol. Chem.* 276: 41594-41602.
6. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 602821. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
7. Macioce, P., et al. 2003.  $\beta$ -dystrobrevin interacts directly with kinesin heavy chain in brain. *J. Cell Sci.* 116: 4847-4856.

## CHROMOSOMAL LOCATION

Genetic locus: KIF5C (human) mapping to 2q23.1.

## SOURCE

NKHC2 (E-5) is a mouse monoclonal antibody raised against amino acids 371-430 mapping within an internal region of NKHC2 of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG<sub>1</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

NKHC2 (E-5) is available conjugated to agarose (sc-374468 AC), 500  $\mu$ g/0.25 ml agarose in 1 ml, for IP; to HRP (sc-374468 HRP), 200  $\mu$ g/ml, for WB, IHC(P) and ELISA; to either phycoerythrin (sc-374468 PE), fluorescein (sc-374468 FITC), Alexa Fluor® 488 (sc-374468 AF488), Alexa Fluor® 546 (sc-374468 AF546), Alexa Fluor® 594 (sc-374468 AF594) or Alexa Fluor® 647 (sc-374468 AF647), 200  $\mu$ g/ml, for WB (RGB), IF, IHC(P) and FCM; and to either Alexa Fluor® 680 (sc-374468 AF680) or Alexa Fluor® 790 (sc-374468 AF790), 200  $\mu$ g/ml, for Near-Infrared (NIR) WB, IF and FCM.

## APPLICATIONS

NKHC2 (E-5) is recommended for detection of NKHC2 of human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2  $\mu$ g per 100-500  $\mu$ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for NKHC2 siRNA (h): sc-61203, NKHC2 shRNA Plasmid (h): sc-61203-SH and NKHC2 shRNA (h) Lentiviral Particles: sc-61203-V.


Molecular Weight of NKHC2: 115 kDa.

Positive Controls: human hippocampus tissue extract or IMR-32 cell lysate: sc-2409.


## RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## DATA



NKHC2 (E-5): sc-374468. Western blot analysis of NKHC2 expression in IMR-32 whole cell lysate.



NKHC2 (E-5): sc-374468. Western blot analysis of NKHC2 expression in human hippocampus tissue extract.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

Alexa Fluor® is a trademark of Molecular Probes, Inc., Oregon, USA