# SANTA CRUZ BIOTECHNOLOGY, INC.

# PBK (C-3): sc-390812



# BACKGROUND

Protein kinases comprise a large group of encoded factors that regulate cellular processes by catalyzing the transfer of a phosphate group to a hydroxyl acceptor in serine, threonine or tyrosine residues. Kinases are capable of influencing the oncogenic potential of cell sytems at the level of oncoprotein or tumor suppressor protein phosphorylation states. Human PDZ-binding kinase, known as PBK, is a 322 amino acid, T/SXV motif-containing serine/threonine kinase that is abuntant in placenta and absent from adult brain tissue. A PDZ domain in the tumor suppressor protein Dlg can coordinate with the T/SXV motif of PBK. The cell cycle checkpoint kinase Cdc2/cyclin B is an upstream effector of PBK that can phosphorylate and activate PBK. Active PBK may associate with PDZ-containing proteins and influence cell cycle control or cellular proliferation.

# REFERENCES

- 1. Hunter, T. 1995. Protein kinases and phosphatases: the yin and yang of protein phosphorylation and signaling. Cell 80: 225-236.
- 2. Hunter, T. 2000. Signaling—2000 and beyond. Cell 100: 113-127.
- Gaudet, S., et al. 2000. Characterization of PDZ-binding kinase, a mitotic kinase. Proc. Natl. Acad. Sci. USA 97: 5167-5172.
- Abe, Y., et al. 2000. Cloning and expression of a novel MAPKK-like protein kinase, lymphokine-activated killer T-cell-originated protein kinase, specifically expressed in the testis and activated lymphoid cells. J. Biol. Chem. 275: 21525-21531.
- Zhao, S., et al. 2001. PDZ-binding kinase participates in spermatogenesis. Int. J. Biochem. Cell Biol. 33: 631-636.
- Dougherty, J.D., et al. 2005. PBK/TOPK, a proliferating neural progenitorspecific mitogen-activated protein kinase kinase. J. Neurosci. 25: 10773-10785.
- Nandi, A.K., et al. 2007. Attenuation of DNA damage checkpoint by PBK, a novel mitotic kinase, involves protein-protein interaction with tumor suppressor p53. Biochem. Biophys. Res. Commun. 358: 181-188.

### **CHROMOSOMAL LOCATION**

Genetic locus: PBK (human) mapping to 8p21.1; Pbk (mouse) mapping to 14 D1.

# SOURCE

PBK (C-3) is a mouse monoclonal antibody raised against amino acids 19-252 mapping near the N-terminus of PBK of human origin.

# PRODUCT

Each vial contains 200  $\mu g$  lgM kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

# **STORAGE**

Store at 4° C, \*\*D0 NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### **APPLICATIONS**

PBK (C-3) is recommended for detection of PBK of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for PBK siRNA (h): sc-106892, PBK siRNA (m): sc-152042, PBK shRNA Plasmid (h): sc-106892-SH, PBK shRNA Plasmid (m): sc-152042-SH, PBK shRNA (h) Lentiviral Particles: sc-106892-V and PBK shRNA (m) Lentiviral Particles: sc-152042-V.

Molecular Weight of PBK: 36 kDa.

Positive Controls: PC-12 cell lysate: sc-2250, HeLa whole cell lysate: sc-2200 or 3T3-L1 cell lysate: sc-2243.

# **RECOMMENDED SUPPORT REAGENTS**

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-lgGκ BP-HRP: sc-516102 or m-lgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein L-Agarose: sc-2336 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-lgGκ BP-FITC: sc-516140 or m-lgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

#### DATA





PBK (C-3): sc-390812. Western blot analysis of PBK expression in HeLa (A), WEHI-231 (B), M1 (C), NIH/3T3 (D), SUP-T1 (E) and Daudi (F) whole cell lysates. PBK (C-3): sc-390812. Western blot analysis of PBK expression in SW480 (A), 3T3-L1 (B), RAW 264.7 (C) and PC-12 (D) whole cell lysates.

#### **RESEARCH USE**

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.