



PETAG, INC.

255 KEYES AVENUE, HAMPSHIRE, IL, 60140

Telephone: 847-683-2288

Fax: 847-683-2003

Customer Service: 800-323-6878

Technical Service: 800-323-0877

Website: www.petag.com

Email: info@petag.com



Every effort has been made to ensure the accuracy of the information published. However, it remains the responsibility of the readers to familiarize themselves with the product information contained on the USA product label or package insert.

GOATS MILK ESBILAC® LIQUID *PetAg*

For Puppies

All Natural Whole Goats Milk, No Preservatives

Milk formula for puppies with sensitive digestive systems

Food supplement for adult dogs

For puppies newborn to 6 weeks of age

All natural with vitamins, minerals and trace nutrients

SHAKE GENTLY

REFRIGERATE AFTER OPENING

Directions for Feeding Goats Milk Esbilac® (GME™): All puppies should receive their mother's milk for at least 2 days, if possible. The colostrum milk gives extra nutrition and temporary immunity against some diseases. Warm GME to room or body temperature. **Feed puppies 2 tablespoons (30 mL) per 4 ounces (115 g) of body weight daily. Divide the daily feeding amount into equal portions for each feeding.** Small and/or weak puppies should be fed every 3 to 4 hours. Larger and/or older puppies can do well being fed every 8 hours. The amount required should be increased or decreased to meet the individual requirements for each puppy. The PetAg 4 oz. Nurser is suited for feeding most puppies. When puppies are old enough to lap, begin offering GME in a shallow container. During the 6th week mix GME with Esbilac® 2nd Step™ Weaning Food to allow the puppy to gradually switch to solid food.

Pregnant and Lactating Bitches: Feed 2 tablespoons (30 mL) GME per 5 lbs. (2.2 kg) body weight daily. **Growing Puppies (post weaning), Show Dogs and/or Convalescing Dogs:** Feed 1 tablespoon (15 mL) GME per 5 lbs. (2.2 kg) body weight daily.

Discard after 72 hours. Do not freeze in original container. Not for human consumption.

Ingredients: Water, whole goats milk, soybean oil, calcium caseinate, calcium carbonate, choline chloride, DL-methionine, L-arginine, lecithin, carrageenan, monopotassium phosphate, potassium chloride, salt, tricalcium phosphate, dicalcium phosphate, magnesium sulfate, zinc sulfate, calcium phosphate, vitamin D3 supplement, copper sulfate, vitamin E supplement, vitamin A supplement, niacin supplement, calcium pantothenate, maltodextrins, thiamine hydrochloride, pyridoxine hydrochloride, riboflavin, manganese sulfate, vitamin B12 supplement, folic acid, menadione sodium bisulfite complex, biotin, calcium iodate.

Guaranteed Analysis

Crude Protein, min.	4.5%
Crude Fat, min.	6.0%
Crude Fiber, max.	0.0%
Moisture, max.	85.0%

Calorie Content (Calculated): The calorie content (ME) is 820 kcal/kg or 12.3 kcal/tbs.

Powdered GME is also available.

©2009 PetAg, Inc., 255 Keyes Avenue, Hampshire, IL 60140 U.S.A.

1-800-323-0877

Visit our website www.petag.com

Made in the U.S.A.

International copyright reserved

Use before date stamped on can

	NET WT.	
12 x 11 FL. OZ. (325 mL)	12 x 11 FL. OZ. (312g)	733 0509

NAC No.: 1097017.8

By using this site, you agree to the [Terms of Service](#) and [Privacy Policy](#).

Copyright 2015 - North American Compendiums