



PETAG, INC.

255 KEYES AVENUE, HAMPSHIRE, IL, 60140

Telephone: 847-683-2288

Fax: 847-683-2003

Customer Service: 800-323-6878

Technical Service: 800-323-0877

Website: www.petag.com

Email: info@petag.com



Every effort has been made to ensure the accuracy of the information published. However, it remains the responsibility of the readers to familiarize themselves with the product information contained on the USA product label or package insert.

GOATS MILK ESBILAC® POWDER

PetAg

For Puppies

All Natural Whole Goats Milk, No Preservatives

Milk formula for puppies with sensitive digestive systems

Food supplement for adult dogs

For puppies newborn to 6 weeks of age

Directions for Mixing Goats Milk Esbilac® (GME™) Powder

Gently stir or shake 1 part powdered GME into 2 parts warm water (a part is a teaspoon, tablespoon or cup). Do not mix more GME than can be consumed in 24 hours. Do not use a blender. Recommended GME should be kept refrigerated. One 12 oz. (340g) can of powder makes approximately 60 oz. (1800mL) of liquid mixed 1:2. One 5 lb. (2.28kg) bag of powder makes approximately 3-1/3 gallons of liquid mixed 1:2.

Directions for Feeding

All puppies should receive their mother's milk for at least 2 days, if possible. This colostrum milk gives extra nutrition and temporary immunity against some diseases. Warm reconstituted GME to room or body temperature before serving. **Feed puppies 2 tablespoons (30mL) of liquid GME per 4 ounces (115g) of body weight daily. The daily feeding amount should be divided into equal portions for each feeding.** Puppies' needs will vary and this amount may have to be increased or decreased depending on the individual. Young and/or weak puppies should be fed every 3 to 4 hours. Older puppies can do well being fed every 8 hours. Weigh the puppies 3 times per week to assure adequate feeding. **Consult your veterinarian for additional advice.** The PetAg 2 oz (60mL) or 4 oz (120mL) Small Animal Nurser is suited for feeding most puppies. When puppies are old enough to lap, begin offering reconstituted GME in a saucer. During the 5th week mix reconstituted GME with Esbilac® 2nd Step™ Puppy Weaning Food. Use of the GME and 2nd Step will allow the puppy to be gradually switched to solid food. **Pregnant and lactating bitches:** Feed 2 teaspoons (4g) GME powder per 5 lbs. (2.2kg) body weight daily. **Growing puppies already on solid food, show dogs, and/or convalescing dogs:** Feed 1 teaspoon (2g) GME powder per 5 lbs. (2.2kg) body weight daily by top dressing food.

Storage

Reconstituted MGE must be kept refrigerated for up to 24 hours. Opened powder must be refrigerated for up to 3 months, or can be frozen for up to 6 months to preserve freshness.

Ingredients

Dried whole goat milk, soybean oil, calcium sodium caseinate, calcium carbonate, choline chloride, monopotassium phosphate, salt, L-arginine, DL-methionine, lecithin, dipotassium phosphate, carrageenan, magnesium sulfate, zinc sulfate, ferrous sulfate, silicon dioxide, copper sulfate, vitamin E supplement, calcium pantothenate, ascorbic acid, niacinamide, maltodextrins, vitamin A supplement, vitamin D3 supplement, manganese sulfate, folic acid, riboflavin, vitamin B12 supplement, thiamine hydrochloride, pyridoxine hydrochloride, biotin.

Guaranteed Analysis

Crude Protein, min	33.0%
Crude Fat, min	40.0%
Crude Fiber, max	0.0%
Moisture, max	5.0%

Calorie Content (Calculated)

The calorie content (ME) is 900 kcal/kg or 13.5 kal/Tbs when powder is reconstituted 1:2.

See back of label for recipes and freezing instructions.

NOT FOR HUMAN CONSUMPTION.

Sold by weight not by volume

All states have laws governing the care of wildlife, and there are a number of potential health hazards that may exist. We have available a list of wildlife professionals nationwide. Call a member of the NWRA, or our toll free number for assistance.

©2008, 2009 PetAg, Inc., 255 Keyes Avenue, Hampshire, IL 60140 U.S.A.

1-800-323-0877

Visit our website: www.petag.com

Made in U.S.A.

International Copyright Reserved.

Use before date stamped on top and bottom. Store in a cool, dry place.

NET WT.	
3/4 oz	
150 g	
12 OZ. (340g)	732 0309
5 lb (2.27 kg)	861F 0308, 861B 0308
12 x emerg. kit: 3/4 oz. powder, .035 oz. Bene-Bac One Shot/case	

NAC No.: 1097038.5

By using this site, you agree to the [Terms of Service](#) and [Privacy Policy](#).

Copyright 2015 - North American Compendiums