

# Mo-CD8/CD38

## 2 Color FCM Reagent: sc-3983



### BACKGROUND

Mouse CD8/CD38: sc-3983 is a direct immunofluorescence reagent formatted to identify and determine the percentage of activated suppressor/cytotoxic T lymphocytes in erythrocyte-lysed whole blood, based on cell-surface antigen expression. CD8 identifies suppressor/cytotoxic T lymphocytes (1,2) and binds class I MHC molecules, resulting in enhanced activation of resting T lymphocytes (3). CD38 is expressed during the early and final stages of T and B cell differentiation, but not during the intermediate stage (4). CD38 is also expressed on activated T and NK lymphocytes (5). CD8+CD38+ T lymphocytes exhibit decreased activity of CD73, which is a maturation marker for T and B lymphocytes, suggesting that these cells are activated and/or immature suppressor/cytotoxic T lymphocytes (6). In HIV infection, the coexpression of CD8 and CD38, with a concomitant low percentage of CD4+ lymphocytes, is closely associated with disease progression (7). Selective elevation of HLA-DR+ CD38-CD8+ cells may be a marker of stable HIV disease (8,9).

Antigen Expression	Cell Type Identified
CD38+	Activated T Lymphocytes
CD8+	Suppressor/Cytotoxic T Cells
CD8+ CD38+	Activated Suppressor/Cytotoxic T Cells

### STORAGE

Store at 4° C. Do not freeze. Stable for one year from the date of shipment. Protect reagents from prolonged exposure to light.

### PRODUCT

Supplied in 1.0 ml of PBS containing 0.1% azide and 0.1% gelatin. Sufficient for 50 tests. This product has been titrated for optimal performance. Recommended use is 20 uL per test (1x10<sup>6</sup> cells). **For research use only. Not for use in diagnostic procedures.**

### INSTRUMENT

Mouse CD8/CD38: sc-3983 is recommended for use with either a single or dual laser Flow Cytometer fitted with appropriate acquisition and analysis software, such as the FACSCalibur™ Flow Cytometer fitted with CellQuest™ Software by Becton Dickinson.

The flow cytometer must be equipped with a 488 nm laser and must be capable of detecting light scatter (forward and side) and two-color fluorescence with emission detectable in two ranges: 515-545 nm, 562-607 nm.

### ISOTYPE CONTROL

sc-3983 CON (rat IgG<sub>2a</sub> FITC/rat IgG<sub>2a</sub> PE) is the isotype matched negative control for this system and is suitable for 50 tests.

Antigen	Clone	Isotype	Label*	Detection Range (nm)
CD8	53-6.7	rat IgG <sub>2a</sub>	FITC	515-545
CD38	90	rat IgG <sub>2a</sub>	PE	562-607

\*Fluorescent labels include FITC: Fluorescein isothiocyanate; PE: phycoerythrin

### REFERENCES

- Ledbetter, J.A., Evans, R.L., Lipinski, M., Cunningham-Rundles, C., Good, R.A., and Herzenberg, L.A. 1981. Evolutionary conservation of surface molecules that distinguish T-lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. *J. Exp. Med.* **153**: 310-323.
- Evans, R.L., Wall, D.W., Platsoucas, C.D., *et al.* 1981. Thymus-dependent membrane antigens in man: Inhibition of cell-mediated lympholysis by monoclonal antibodies to the TH2 antigen. *Proc. Natl. Acad. Sci. USA* **78**: 544-548.
- Gallagher, P.F., Fazekas de St. Groth, B., and Miller, J.F. 1989. CD4 and CD8 molecules can physically associate with the same T-cell receptor. *Proc. Natl. Acad. Sci. USA* **86**: 10044-10048.
- Dörken, B., Möller, P., Pezzutto, A., Schwartz-Albiez, R., and Moldenhauer, G. B-cell antigens: CD19. In: Knapp, W., Dörken, B., Gilks, W.R. *et al.* eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989: 34-36.
- Sconocchia, G., Titus, J.A., Mazzoni, A., Visintin, A., Pericle, F., Hicks, S.W., Malavasi, F., and Segal, D.M. 1999. CD38 triggers cytotoxic responses in activated human natural killer cells. *Blood* **94**: 3864-3871.
- Salazar-Gonzalez, J.F., Moody, D.J., Giorgi, J.V., Martinez-Maza, O., Mitsuyasu, R.T., and Fahey, J.L. 1985. Reduced ecto-5'-nucleotidase activity and enhanced OKT10 and HLA-DR expression of CD8 (T suppressor/cytotoxic) lymphocytes in the acquired immune deficiency syndrome: Evidence of CD8 cell immaturity. *J. Immunol.* **135**: 1778-1785.
- Giorgi, J. and Hultin, L. 1990. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. *Clin. Immunol. Newslett.* **10**: 55-62.
- Liu, Z., Cumberland, W.G., Hultin, L.E., Prince, H.E., Detels, R., and Giorgi, J.V. 1997. Elevated CD38 antigen expression on CD8+ T cells is a stronger marker for the risk of chronic HIV disease progression to AIDS and death in the Multicenter AIDS Cohort Study than CD4+ cell count, soluble immune activation markers, or combinations of HLA-DR and CD38 expression. *J. Acquir. Immune Defic. Syndr. Hum. Retrovirol.* **16**: 83-92.
- Giorgi, J.V., Ho, H.N., Hirji, K., Chou, C.C., *et al.* 1994. CD8+ lymphocyte activation at human immunodeficiency virus type 1 seroconversion: development of HLA-DR+ CD38- CD8+ cells is associated with subsequent stable CD4+ cell levels. The Multicenter AIDS Cohort Study Group. *J. Infect. Dis.* **170**: 775-781.