

# Cdk7 (1-346): sc-4103 WB

## BACKGROUND

Progression through the cell cycle requires activation of a series of enzymes designated cyclin dependent kinases (Cdks). The monomeric catalytic subunit Cdk2, a critical enzyme for initiation of cell cycle progression, is completely inactive. Partial activation is achieved by the binding of regulatory cyclins such as cyclin D1, while full activation requires additional phosphorylation at Thr 160. The enzyme responsible for the phosphorylation of Cdk2 on Thr 160 and also of Cdc2 p34 on Thr 161, designated Cdk-activating kinase (CAK), has been partially purified and shown to be comprised of a 42 kDa catalytic subunit and a 37 kDa regulatory subunit. The catalytic subunit, designated Cdk7, has been identified as the mammalian homolog of MO15, a protein kinase demonstrated in starfish and *Xenopus*. The regulatory subunit is a novel cyclin (cyclin H) and is required for activation of Cdk7. Like other Cdks, Cdk7 contains a conserved threonine residue required for full activity; mutation of this residue severely reduces CAK activity.

## REFERENCES

1. Nurse, P. 1994. Ordering S phase and M phase in the cell cycle. *Cell* 79: 547-550.
2. Sherr, C.J. 1994. G1 phase progression: cycling on cue. *Cell* 79: 551-555.
3. Helchman, K.A. and Roberts, J.M. 1994. Rules to replicate by. *Cell* 79: 557-562.
4. King, R.W., Jackson, P.K., and Kirschner, M.W. 1994. Mitosis in transition. *Cell* 79: 563-571.
5. Hunter, T. and Pines, J. 1994. Cyclins and cancer II: cyclin D and Cdk inhibitors come of age. *Cell* 79: 573-582.
6. Levedakou, E.N., He, M., Baptist, E.W., Craven, R.J., Cance, W.G., Welch, P.L., Simmons, A., Naylor, S.L., Leach, R.J., Lewis, T.B., Bowcock, A., Liu, E.T. 1994. Two novel human serine/threonine kinases with homologies to the cell cycle regulating *Xenopus* MO15, and NIMA kinases: cloning and characterization of their expression pattern. *Oncogene* 9: 1977-1988.
7. Wu, L., Yee, A., Liu, L., Carbonaro-Hall, D., Venkatesan, N., Tolo, V.T., and Hall, F.L. 1994. Molecular cloning of the human CAK1 gene encoding a cyclin-dependent kinase-activating kinase. *Oncogene* 9: 2089-2096.
8. Kato, J.Y., Matsuoka, M., Strom, D.K., and Sherr, C.J. 1994. Regulation of cyclin D-dependent kinase 4 (Cdk4) by Cdk4-activating kinase. *Mol. Cell Biol.* 14: 2713-2721.
9. Matsuoka, M., Kato, J.Y., Fisher, R.P., Morgan, D.O., and Sherr, C.J. 1994. Activation of cyclin-dependent kinase 4 (Cdk4) by mouse MO15-associated kinase. *Mol. Cell Biol.* 14: 7265-7275.

## CHROMOSOMAL LOCATION

Genetic locus: CDK7 (human) mapping to 5q12.1; Cdk7 (mouse) mapping to 13 D1.

## SOURCE

Cdk7 (1-346) is expressed in *E. coli* as a 42 kDa polyhistidine tagged fusion protein corresponding to amino acids 1-346 representing full length Cdk7 of human origin.

## PRODUCT

Cdk7 (1-346) is purified from bacterial lysates (> 98%) by Ni<sup>++</sup> affinity chromatography; supplied as 10 µg in 0.1 ml SDS-PAGE loading buffer.

## APPLICATIONS

Cdk7 (1-346) is suitable as a Western blotting control for sc-529, sc-723, sc-856 and sc-7344.

## SELECT PRODUCT CITATIONS

1. Acevedo-Duncan, M., Patel, R., Whelan, S. and Bicaku E. 2002. Human glioma PKC- $\alpha$  and PKC- $\beta$ II phosphorylate cyclin-dependent kinase activating kinase during the cell cycle. *Cell Prolif.* 35: 23-36.

## STORAGE

Store -20° C; stable for one year from the date of shipment.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.