

# APG5 siRNA (m): sc-41446

## BACKGROUND

In yeast, autophagy is an essential process for survival during nutrient starvation and cell differentiation. The process of autophagy is characterized as a non-selective degradation of cytoplasmic proteins into membrane structures called autophagosomes, and it is dependent on several proteins, including the autophagy proteins ATG5 (APG5) and APG7. Yeast Apg7 and the human homolog, APG7, share similarities with the ubiquitin-activating enzyme E1 in *Saccharomyces cerevisiae* and are likewise responsible for enzymatically activating the autophagy conjugation system. APG5 and the human homolog, ATG5, also designated APG5, apoptosis-specific protein or APS, function as substrates for the autophagy protein Apg12. These proteins are covalently bonded together to form Apg12/ATG5 conjugates, which are required for the progression of autophagy.

## CHROMOSOMAL LOCATION

Genetic locus: Atg5 (mouse) mapping to 10 B2.

## PRODUCT

APG5 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see APG5 shRNA Plasmid (m): sc-41446-SH and APG5 shRNA (m) Lentiviral Particles: sc-41446-V as alternate gene silencing products.

For independent verification of APG5 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-41446A, sc-41446B and sc-41446C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

APG5 siRNA (m) is recommended for the inhibition of APG5 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

APG5 (C-1): sc-133158 is recommended as a control antibody for monitoring of APG5 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor APG5 gene expression knockdown using RT-PCR Primer: APG5 (m)-PR: sc-41446-PR (20  $\mu$ l, 449 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## SELECT PRODUCT CITATIONS

- Yang, C., et al. 2008. Autophagy is associated with apoptosis in cisplatin injury to renal tubular epithelial cells. *Am. J. Physiol. Renal Physiol.* 294: F777-F787.
- Ahn, J.H., et al. 2011. Suppression of autophagy sensitizes multidrug resistant cells towards Src tyrosine kinase specific inhibitor PP2. *Cancer Lett.* 310: 188-197.
- Chung, Y.H., et al. 2012. Microtubule-associated protein light chain 3 regulates Cdc42-dependent Actin ring formation in osteoclast. *Int. J. Biochem. Cell Biol.* 44: 989-997.
- Liao, C.C., et al. 2013. Recombinant protein rVP1 upregulates BECN1-independent autophagy, MAPK1/3 phosphorylation and MMP9 activity via WIPI1/WIPI2 to promote macrophage migration. *Autophagy* 9: 5-19.
- Yang, Z., et al. 2014. Hypoxia induces autophagic cell death through hypoxia-inducible factor 1 $\alpha$  in microglia. *PLoS ONE* 9: e96509.
- Ozeki, N., et al. 2015. Interleukin-1 $\beta$ -induced autophagy-related gene 5 regulates proliferation of embryonic stem cell-derived odontoblastic cells. *PLoS ONE* 10: e0124542.
- Sin, J., et al. 2016. Mitophagy is required for mitochondrial biogenesis and myogenic differentiation of C2C12 myoblasts. *Autophagy* 12: 369-380.
- Ye, J., et al. 2017. The role of autophagy in pro-inflammatory responses of microglia activation via mitochondrial reactive oxygen species *in vitro*. *J. Neurochem.* 142: 215-230.
- Xiao, Y., et al. 2018. TAX1BP1 overexpression attenuates cardiac dysfunction and remodeling in STZ-induced diabetic cardiomyopathy in mice by regulating autophagy. *Biochim. Biophys. Acta Mol. Basis Dis.* 1864: 1728-1743.
- Seong, S.B., et al. 2019. Autophagy precedes apoptosis in Angiotensin II-induced podocyte injury. *Cell. Physiol. Biochem.* 53: 747-759.
- Sivangala Thandi, R., et al. 2020. Ornithine-A urea cycle metabolite enhances autophagy and controls *Mycobacterium tuberculosis* infection. *Nat. Commun.* 11: 3535.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.