

P2X7 siRNA (h): sc-42575

BACKGROUND

The P2X receptor family is comprised of ligand-gated ion channels that allow for the increased permeability of calcium into the cell in response to extracellular ATP. The seven P2X receptors, P2X1-P2X7, form either homomeric or heteromeric channels or both. They are characterized by intracellular amino- and carboxy-termini. P2X receptors are expressed in a wide variety of tissues, including neurons, prostate, bladder, pancreas, colon, testis and ovary. The major function of the P2X receptors is to mediate synaptic transmissions between neurons and to other tissues via the binding of extracellular ATP, which acts as a neurotransmitter. The P2X receptors may be involved in the onset of necrosis or apoptosis after prolonged exposure to high concentrations of extracellular ATP.

CHROMOSOMAL LOCATION

Genetic locus: P2RX7 (human) mapping to 12q24.31.

PRODUCT

P2X7 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see P2X7 shRNA Plasmid (h): sc-42575-SH and P2X7 shRNA (h) Lentiviral Particles: sc-42575-V as alternate gene silencing products.

For independent verification of P2X7 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-42575A, sc-42575B and sc-42575C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

P2X7 siRNA (h) is recommended for the inhibition of P2X7 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

P2X7 (D-1): sc-514962 is recommended as a control antibody for monitoring of P2X7 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor P2X7 gene expression knockdown using RT-PCR Primer: P2X7 (h)-PR: sc-42575-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

- Deli, T., et al. 2007. Functional genomics of calcium channels in human melanoma cells. *Int. J. Cancer* 121: 55-65.
- Chen, J.B., et al. 2012. Adenosine-5'-triphosphate up-regulates proliferation of human cardiac fibroblasts. *Br. J. Pharmacol.* 166: 1140-1150.
- Hwang, Y.S., et al. 2014. Lysophosphatidic acid stimulates osteoclast fusion through OC-STAMP and P2X7 receptor signaling. *J. Bone Miner. Metab.* 32: 110-122.
- Ghalali, A., et al. 2017. Atorvastatin decreases HBx-induced phospho-Akt in hepatocytes via P2X receptors. *Mol. Cancer Res.* 15: 714-722.
- Zheng, H., et al. 2017. ATM-activated autotaxin (ATX) propagates inflammation and DNA damage in lung epithelial cells: a new mode of action for silica-induced DNA damage? *Carcinogenesis* 38: 1196-1206.
- Liu, Y.H., et al. 2018. The ATP-P2X7 signaling axis is an essential sentinel for intracellular *Clostridium difficile* pathogen-induced inflammasome activation. *Front. Cell. Infect. Microbiol.* 8: 84.
- Mawatwal, S., et al. 2018. Calcimycin induced IL-12 production inhibits intracellular mycobacterial growth by enhancing autophagy. *Cytokine* 111: 1-12.
- Cao, Y., et al. 2019. Extracellular and macropinocytosis internalized ATP work together to induce epithelial-mesenchymal transition and other early metastatic activities in lung cancer. *Cancer Cell Int.* 19: 254.
- Ghalali, A., et al. 2020. PTEN and PHLPP crosstalk in cancer cells and in TGF β -activated stem cells. *Biomed. Pharmacother.* 127: 110112.
- Mishra, A., et al. 2021. Soybean lectin induces autophagy through P2RX7 dependent activation of NF κ B-ROS pathway to kill intracellular mycobacteria. *Biochim. Biophys. Acta Gen. Subj.* 1865: 129806.
- Zhong, W., et al. 2022. Nanocarrier-assisted delivery of metformin boosts remodeling of diabetic periodontal tissue via cellular exocytosis-mediated regulation of endoplasmic reticulum homeostasis. *ACS Nano* 16: 19096-19113.

RESEARCH USE

For research use only, not for use in diagnostic procedures.