

CHIP siRNA (h): sc-43555

BACKGROUND

CHIP (carboxy terminus of HSP 70-interacting protein), also designated STIP1 homology and U-box containing protein 1, HSPABP2, NY-CO-7, SDCCAG7 and STUB1, is a cytoplasmic E3 ubiquitin ligase that influences protein ubiquitylation. CHIP interacts with Smad1/Smad4 and blocks BMP signaling through the ubiquitin-mediated degradation of Smad proteins. CHIP controls both association of HSP70/Hsp90 chaperones with ErbB2 and down-regulation of ErbB2 induced by inhibitors of Hsp90. A 1.3-kb transcript is most abundant in striated muscle (heart and skeletal muscle), with lower expression in pancreas and brain.

CHROMOSOMAL LOCATION

Genetic locus: STUB1 (human) mapping to 16p13.3.

PRODUCT

CHIP siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CHIP shRNA Plasmid (h): sc-43555-SH and CHIP shRNA (h) Lentiviral Particles: sc-43555-V as alternate gene silencing products.

For independent verification of CHIP (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-43555A, sc-43555B and sc-43555C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

CHIP siRNA (h) is recommended for the inhibition of CHIP expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

CHIP (G-2): sc-133066 is recommended as a control antibody for monitoring of CHIP gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CHIP gene expression knockdown using RT-PCR Primer: CHIP (h)-PR: sc-43555-PR (20 μ l, 416 bp). Annealing temperature for the primers should be $55-60^{\circ}$ C and the extension temperature should be $68-72^{\circ}$ C.

SELECT PRODUCT CITATIONS

- Kumar, P., et al. 2007. CHIP and HSPs interact with β -APP in a proteasome-dependent manner and influence A β metabolism. *Hum. Mol. Genet.* 16: 848-864.
- Jeong, J.H., et al. 2008. Quercetin-induced ubiquitination and down-regulation of Her-2/neu. *J. Cell. Biochem.* 105: 585-595.
- Wang, J., et al. 2011. Gambogic acid-induced degradation of mutant p53 is mediated by proteasome and related to CHIP. *J. Cell. Biochem.* 112: 509-519.
- Blessing, N.A., et al. 2014. The E3 ligase CHIP mediates ubiquitination and degradation of mixed-lineage kinase 3. *Mol. Cell. Biol.* 34: 3132-3143.
- Zhang, N., et al. 2016. Glycogen synthase kinase-3 β antagonizes Ros-inducible hepatocellular carcinoma cell death through suppression of the apoptosis signal-regulating kinase 1. *Med. Oncol.* 33: 60.
- Wang, J., et al. 2016. PC-1 works in conjunction with E3 ligase CHIP to regulate androgen receptor stability and activity. *Oncotarget* 7: 81377-81388.
- Shao, M., et al. 2016. E3 ubiquitin ligase CHIP interacts with C-type lectin-like receptor CLEC-2 and promotes its ubiquitin-proteasome degradation. *Cell. Signal.* 28: 1530-1536.
- Yong, H.J., et al. 2017. Von Hippel-Lindau regulates interleukin-32 β stability in ovarian cancer cells. *Oncotarget* 8: 69833-69846.
- Cho, S.H., et al. 2017. The antitumor effect of C-terminus of HSP 70-interacting protein via degradation of c-Met in small cell lung cancer. *Korean J. Thorac. Cardiovasc. Surg.* 50: 153-162.
- Blessing, N.A., et al. 2017. Osmotic and heat stress-dependent regulation of MLK4 β and MLK3 by the CHIP E3 ligase in ovarian cancer cells. *Cell. Signal.* 39: 66-73.
- Sun, L.W., et al. 2021. Corosolic acid attenuates the invasiveness of glioblastoma cells by promoting CHIP-mediated Axl degradation and inhibiting GAS6/Axl/JAK Axis. *Cells* 10: 2919.

RESEARCH USE

For research use only, not for use in diagnostic procedures.