

mAChR M5 (230-426): sc-4506 WB

BACKGROUND

The muscarinic acetylcholine receptors (mAChR) mediate a variety of cellular responses, including inhibition of adenylate cyclase, breakdown of phosphoinositides and modulation of potassium channels. The mAChRs transduce signals by coupling to G proteins, which then modulate several downstream effector proteins and ion channels. Five mAChR subtypes have been identified, designated M1 to M5. The five receptor subtypes show distinct patterns of tissue distribution, as well as distinct pharmacological and functional properties. The amino acid sequence of each mAChR subtype reflects a structure that is characteristic of G protein-coupled receptors, consisting of seven highly conserved transmembrane segments and a large intracellular region unique to each subtype, which constitutes the effector-coupling domain

REFERENCES

1. Peralta, E.G., Winslow, J.W., Peterson, G.L., Smith, D.H., Ashkenazi, A., Ramachandran, J., Schimerlik, M.I., and Capon, D.J. 1987. Primary structure and biochemical properties of an M2 muscarinic receptor. *Science* 236: 600-605.
2. Liao, C.F., Themmen, A.P., Joho, R., Barberis, C., Birnbaumer, M., and Birnbaumer, L. 1989. Molecular cloning and expression of a fifth muscarinic acetylcholine receptor. *J. Biol. Chem.* 264: 7328-7337.
3. Hulme, E.C. 1990. Muscarinic acetylcholine receptors: typical G-coupled receptors. *Symp. Soc. Exp. Biol.* 44: 39-54.
4. Hulme, E.C., Kurtenbach, E., and Curtis, C.A. 1991. Muscarinic acetylcholine receptors: structure and function. *Biochem. Soc. Trans.* 19: 133-138.
5. Caulfield, M.P. 1993. Muscarinic receptor-characterization, coupling and function. *Pharmacol. Ther.* 58: 319-379.
6. Brann, M.R., Ellis, J., Jorgensen, H., Hill-Eubanks, D., and Jones, S.V. 1993. Muscarinic acetylcholine receptor subtypes: localization and structure/function. *Prog. Brain Res.* 98: 121-127.
7. Brauner-Osborne, H. and Brann, M.R. 1996. Pharmacology of muscarinic acetylcholine receptor subtypes (M1–M5): high throughput assays in mammalian cells. *Eur. J. Pharmacol.* 295: 93-102.
8. Tice, M.A., Hashemi, T., Taylor, L.A., and McQuade, R.D. 1996. Distribution of muscarinic receptor subtypes in rat brain from postnatal to old age. *Brain Res. Dev. Brain Res.* 92: 70-76.
9. Madziva, M.T., Bai, J., Bhalla, A., Chapman, E.R., and Edwardson, J.M. 2005. Effects of synaptotagmin reveal two distinct mechanisms of agonist-stimulated internalization of the M4 muscarinic acetylcholine receptor. *Br. J. Pharmacol.* 144: 761-771.

SOURCE

mAChR M5 (230-426) is expressed in *E. coli* as a 49 kDa tagged fusion protein corresponding to amino acids 230-426 of mAChR M5 of human origin.

STORAGE

Store at -20° C; stable for one year from the date of shipment.

PRODUCT

mAChR M5 (230-426) is purified from bacterial lysates (>98%) by glutathione agarose affinity chromatography; supplied as 10 µg in 0.1 ml SDS-PAGE loading buffer.

APPLICATIONS

mAChR M5 (230-426) is suitable as a Western blotting control for sc-9110.

RESEARCH USE

For research use only, not for use in diagnostic procedures.