

RPA 14 kDa subunit (M-18): sc-46508

BACKGROUND

The single-stranded-DNA-binding proteins (SSBs) are essential for DNA function in prokaryotic and eukaryotic cells, mitochondria, phages and viruses. Replication protein A (RPA), a highly conserved eukaryotic protein, is a heterotrimeric SSB that is composed of three subunits, designated RPA 14 kDa (also known as RPA3), RPA 32 kDa and RPA 70 kDa. Together, these subunits play an important role in DNA replication, recombination and repair. RPA is one of the major damage-recognition structures involved in the early stage of nucleotide excision repair and may play a role in telomere maintenance. The binding of human RPA (hRPA) to DNA involves molecular polarity, in which initial hRPA binding occurs on the 5' side of a ssDNA substrate and then extends in the 3' direction to create a stably bound hRPA. The RPA 14 kDa subunit localizes to the nucleus and is the smallest component of the RPA complex, functioning with the other subunits to regulate various aspects of DNA metabolism.

CHROMOSOMAL LOCATION

Genetic locus: Rpa3 (mouse) mapping to 6 A1.

SOURCE

RPA 14 kDa subunit (M-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of RPA 14 kDa subunit of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-46508 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-46508 X, 200 µg/0.1 ml.

APPLICATIONS

RPA 14 kDa subunit (M-18) is recommended for detection of RPA 14 kDa subunit of mouse origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for RPA 14 kDa subunit siRNA (m): sc-45713, RPA 14 kDa subunit shRNA Plasmid (m): sc-45713-SH and RPA 14 kDa subunit shRNA (m) Lentiviral Particles: sc-45713-V.

RPA 14 kDa subunit (M-18) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.


Molecular Weight of RPA: 14 kDa.

Positive Controls: NIH/3T3 nuclear extract: sc-2138 or Sol8 nuclear extract: sc-2157.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

DATA



RPA 14 kDa subunit (M-18): sc-46508. Western blot analysis of RPA 14 kDa subunit expression in NIH/3T3 (A), Sol8 (B), HL-60 (C) and K-562 (D) nuclear extracts. Note lack of reactivity in lanes C and D.

SELECT PRODUCT CITATIONS

1. McJunkin, K., et al. 2011. Reversible suppression of an essential gene in adult mice using transgenic RNA interference. Proc. Natl. Acad. Sci. USA 108: 7113-7118.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **RPA 14 kDa subunit (A-2): sc-393891**, our highly recommended monoclonal alternative to RPA 14 kDa subunit (M-18).