

mGluR-3 (N-12): sc-47139

BACKGROUND

The mGluR proteins (metabotropic glutamate receptors) are members of the G protein-coupled receptor family and are functionally and pharmacologically distinct from the GluR proteins (ionotropic glutamate receptors). The eight currently known mGluR proteins are mediated by two G proteins with opposing regulation of adenylate cyclase pathways. The activities of mGluR-1 and mGluR-5 are mediated by a G protein that activates a phosphatidylinositol-calcium second messenger system and generates a calcium-activated chloride current. The remainder of the eight subtypes of mGluR have an activity mediated by a G protein that inhibits adenylate cyclase activity. mGluR-3, which may interact with GRASP, acts as a receptor for glutamate.

REFERENCES

1. Makoff, A., Volpe, F., Lelchuk, R., Harrington, K. and Emson, P. 1997. Molecular characterization and localization of human metabotropic glutamate receptor type 3. *Brain Res. Mol. Brain Res.* 40: 55-63.
2. Kammermeier, P.J. and Yun, J. 2005. Activation of metabotropic glutamate receptor 1 dimers requires glutamate binding in both subunits. *J. Pharmacol. Exp. Ther.* 312: 502-508.
3. Bäckström, P. and Hyttiä, P. 2005. Suppression of alcohol self-administration and cue-induced reinstatement of alcohol seeking by the mGlu2/3 receptor agonist LY379268 and the mGlu8 receptor agonist (S)-3,4-DCPG. *Eur. J. Pharmacol.* 528: 110-108.
4. Pacheco Ojalora, L.F., Cough, J., Shigamoto, R., Zarei, M.M. and Garrido Sanabria, E.R. 2006. Abnormal mGluR-2/-3 expression in the perforant path termination zones and mossy fibers of chronically epileptic rats. *Brain Res.* 1098: 170-185.
5. Yoshimizu, T., Shimazaki, T., Ito, A. and Chaki, S. 2006. An mGluR-2/-3 antagonist, MGS0039, exerts antidepressant and anxiolytic effects in behavioral models in rats. *Psychopharmacology* 186: 587-593.
6. Ohana, L., Barchad, O., Parnas, I. and Parnas, H. 2006. The metabotropic glutamate G protein-coupled receptors mGluR-3 and mGluR-1a are voltage sensitive. *J. Biol. Chem.* 281: 24204-24215.
7. Marengo, S., Steele, S.U., Egan, M.F., Goldberg, T.E., Straub, R.E., Sharrief, A.Z. and Weinberger, D.R. 2006. Effect of metabotropic glutamate receptor 3 genotype on N-acetylaspartate measures in the dorsolateral prefrontal cortex. *Am. J. Psychiatry* 163: 740-742.
8. Bonsi, P., Cuomo, D., Picconi, B., Sciamanna, G., Tschertner, A., Tolu, M., Bernardi, G., Calabresi, P. and Pisani, A. 2006. Striatal metabotropic glutamate receptors as a target for pharmacotherapy in Parkinson's disease. *Amino Acids* 32: 189-95.

CHROMOSOMAL LOCATION

Genetic locus: GRM3 (human) mapping to 7q21.11; Grm3 (mouse) mapping to 5 A1.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

SOURCE

mGluR-3 (N-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an N-terminal extracellular domain of mGluR-3 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-47139 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

mGluR-3 (N-12) is recommended for detection of mGluR-3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

mGluR-3 (N-12) is also recommended for detection of mGluR-3 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for mGluR-3 siRNA (h): sc-61030, mGluR-3 siRNA (m): sc-61031, mGluR-3 shRNA Plasmid (h): sc-61030-SH, mGluR-3 shRNA Plasmid (m): sc-61031-SH, mGluR-3 shRNA (h) Lentiviral Particles: sc-61030-V and mGluR-3 shRNA (m) Lentiviral Particles: sc-61031-V.

Molecular Weight of mGluR-3: 110/220 kDa.

Positive Controls: SK-N-MC cell lysate: sc-2237.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **mGluR-3 (A-10): sc-271899**, our highly recommended monoclonal alternative to mGluR-3 (N-12).