

mGluR-6 (C-16): sc-47148

BACKGROUND

The mGluR proteins (metabotropic glutamate receptors) are members of the G protein-coupled receptor family and are functionally and pharmacologically distinct from the GluR proteins (ionotropic glutamate receptors). The eight currently known mGluR proteins are mediated by two G proteins with opposing regulation of adenylate cyclase pathways. The activities of mGluR-1 and mGluR-5 are mediated by a G protein that activates a phosphatidylinositol-calcium second messenger system and generates a calcium-activated chloride current. The remainder of the eight subtypes of mGluR have an activity mediated by a G protein that inhibits adenylate cyclase activity. mGluR-6 is expressed in the synapses of bipolar cell dendrites. This receptor is involved in mediating synaptic transmission from rod and cone photoreceptors to other neurons.

REFERENCES

1. Hashimoto, T., Inazawa, J., Okamoto, N., Tagawa, Y., Bessho, Y., Honda, Y. and Nakanishi, S. 1997. The whole nucleotide sequence and chromosomal localization of the gene for human metabotropic glutamate receptor subtype 6. *Eur. J. Neurosci.* 9: 1226-1235.
2. Dryja, T.P., McGee, T.L., Berson, E.L., Fishman, G.A., Sandberg, M.A., Alexander, K.R., Derlacki, D.J. and Rajagopalan, A.S. 2005. Night blindness and abnormal cone electroretinogram ON responses in patients with mutations in the GRM6 gene encoding mGluR-6. *Proc. Natl. Acad. Sci. USA* 102: 4884-4889.
3. Zeitz, C., van Genderen, M., Neidhardt, J., Luhmann, U.F., Hoeben, F., Forster, U., Wycisk, K., Mátyás, G., Hoyng, C.B., Riemsdag, F., Meire, F., Cremers, F.P. and Berger, W. 2005. Mutations in GRM6 cause autosomal recessive congenital stationary night blindness with a distinctive scotopic 15-Hz flicker electroretinogram. *Invest. Ophthalmol. Vis. Sci.* 46: 4328-4335.
4. Yang, Z.Q. 2005. Agonists and antagonists for group III metabotropic glutamate receptors 6, 7 and 8. *Curr. Top. Med. Chem.* 5: 913-918.
5. O'Connor, E., Allen, L.E., Bradshaw, K., Boylan, J., Moore, A.T. and Trump, D. 2006. Congenital stationary night blindness associated with mutations in GRM6 encoding glutamate receptor mGluR-6. *Br. J. Ophthalmol.* 90: 653-654.

CHROMOSOMAL LOCATION

Genetic locus: GRM6 (human) mapping to 5q35.3.

SOURCE

mGluR-6 (C-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a C-terminal cytoplasmic domain of mGluR-6 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-47148 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

mGluR-6 (C-16) is recommended for detection of mGluR-6 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

mGluR-6 (C-16) is also recommended for detection of mGluR-6 in additional species, including equine, canine and bovine.

Suitable for use as control antibody for mGluR-6 siRNA (h): sc-61036, mGluR-6 shRNA Plasmid (h): sc-61036-SH and mGluR-6 shRNA (h) Lentiviral Particles: sc-61036-V.

Molecular Weight of mGluR-6: 190 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
 Satisfaction
 Guaranteed

Try **mGluR-6 (1A11): sc-517076**, our highly recommended monoclonal alternative to mGluR-6 (C-16).