

Jumonji (C-17): sc-47604

BACKGROUND

Jumonji, a nuclear protein crucial for neural tube formation, plays an important role in down-regulating atrial natriuretic factor (ANF) gene expression through its interaction with GATA-4 and NKX2-5. Required for normal development of the heart, Jumonji participates in negatively regulating the signaling involved in cell proliferation. Jumonji, which also is designated ARID domain-containing protein 2, is a member of the JMJ transcription factor family of proteins. During embryogenesis Jumonji, which localizes to the nucleus, is expressed primarily in neurons, especially dorsal root ganglion cells.

REFERENCES

1. Berge-Lefranc, J.L., Jay, P., Massacrier, A., Cau, P., Mattei, M.G., Bauer, S., Marsollier, C., Berta, P. and Fontes, M. 1996. Characterization of the human Jumonji gene. *Hum. Mol. Genet.* 5: 1637-1641.
2. Kim, T.G., Chen, J., Sadoshima, J. and Lee, Y. 2004. Jumonji represses atrial natriuretic factor gene expression by inhibiting transcriptional activities of cardiac transcription factors. *Mol. Cell. Biol.* 24: 10151-10160.
3. Jung, J., Kim, T.G., Lyons, G., Kim, H.R. and Lee, Y. 2005. Jumonji regulates cardiomyocyte proliferation via interaction with retinoblastoma protein. *J. Biol. Chem.* 280: 30916-30923.
4. Kim, T.G., Jung, J., Mysliwiec, M.R., Kang, S. and Lee, Y. 2005. Jumonji represses α -cardiac myosin heavy chain expression via inhibiting MEF2 activity. *Biochem. Biophys. Res. Commun.* 329: 544-553.
5. Jensen, L.R., Amende, M., Gurok, U., Moser, B., Gimmel, V., Tzschach, A., Janecke, A.R., Tariverdian, G., Chelly, J., Fryns, J.P., Van Esch, H., Hamel, B., Kleefstra, T., Moraine, C., Gecz, J., Turner, G., Reinhardt, R., Ropers, H.H., Kalscheuer, V.M. and Lenzner, S. 2005. Mutations in the JARID1C gene, which is involved in transcriptional regulation and chromatin remodeling, cause X-linked mental retardation. *Am. J. Hum. Genet.* 76: 227-236.
6. Jung, J., Mysliwiec, M.R. and Lee, Y. 2005. Roles of Jumonji in mouse embryonic development. *Dev. Dyn.* 232: 21-32.

CHROMOSOMAL LOCATION

Genetic locus: JARID2 (human) mapping to 6p22.3; Jarid2 (mouse) mapping to 13 A5.

SOURCE

Jumonji (C-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of Jumonji of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-47604 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-47604 X, 200 μ g/0.1 ml.

APPLICATIONS

Jumonji (C-17) is recommended for detection of Jumonji of mouse and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Jumonji (C-17) is also recommended for detection of Jumonji in additional species, including equine.

Suitable for use as control antibody for Jumonji siRNA (h): sc-60872, Jumonji siRNA (m): sc-60873, Jumonji shRNA Plasmid (h): sc-60872-SH, Jumonji shRNA Plasmid (m): sc-60873-SH, Jumonji shRNA (h) Lentiviral Particles: sc-60872-V and Jumonji shRNA (m) Lentiviral Particles: sc-60873-V.

Jumonji (C-17) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of Jumonji: 150-160 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.