

PARN (N-20): sc-47619

BACKGROUND

Exonucleolytic degradation of the poly(A) tail often initiates the first step in the decay of eukaryotic mRNAs. Poly(A)-specific ribonuclease (PARN), a highly poly(A)-specific 3'-exoribonuclease, efficiently degrades mRNA poly(A) tails. PARN, which also may be designated deadenylating nuclease, may also be involved in nonsense-mediated mRNA decay, a critical process of selective degradation of mRNAs that contain premature stop codons, and in the degradation of inherently unstable mRNAs that contain au-rich elements (AREs) in their 3' untranslated regions. PARN, which can form a homodimer, interacts with KHSRP and can be found in a mRNA decay complex with RENT1, RENT2 and RENT3B. It localizes mainly to the nucleus (may be detected in the nucleolus), but may also localize to the cytoplasm.

REFERENCES

1. Buiting, K., Körner, C., Ulrich, B., Wahle, E. and Horsthemke, B. 2000. The human gene has a truncated copy in the Prader-Willi/Angelman syndrome region on 15q11-q13. *Cytogenet. Cell Genet.* 87: 125-131.
2. Martinez, J., Virtanen, A., Thuresson, A.C., Hellman, U., Astrom, J. and Ren, Y.G. 2000. A 54 kDa fragment of the PARN is an oligomeric, processive, and cap-interacting Poly(A)-specific 3' exonuclease. *J. Biol. Chem.* 275: 24222-24230.
3. Scherl, A., Coute, Y., Deon, C., Calle, A., Hochstrasser, D., Kindbeiter, K., Sanchez, J.C., Greco, A. and Diaz, J.J. 2002. Functional proteomic analysis of human nucleolus. *Mol. Biol. Cell* 13: 4100-4109.
4. Ren, Y.G., Kirsebom, L.A. and Virtanen, A. 2004. Coordination of divalent metal ions in the active site of PARN. *J. Biol. Chem.* 279: 48702-48706.
5. Seal, R., Temperley, R., Wilusz, J., Lightowlers, R.N. and Chrzanowska-Lightowlers, Z.M. 2005. Serum-deprivation stimulates cap-binding by PARN at the expense of eIF4E, consistent with the observed decrease in mRNA stability. *Nucleic Acids Res.* 33: 376-387.
6. Wu, M., Reuter, M., Lilie, H., Liu, Y., Wahle, E. and Song, H. 2005. Structural insight into poly(A) binding and catalytic mechanism of human PARN. *EMBO J.* 24: 4082-4093.

CHROMOSOMAL LOCATION

Genetic locus: PARN (human) mapping to 16p13.12; Parn (mouse) mapping to 16 A1.

SOURCE

PARN (N-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of PARN of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-47619 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-47619 X, 200 µg/0.1 ml.

APPLICATIONS

PARN (N-20) is recommended for detection of PARN of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PARN (N-20) is also recommended for detection of PARN in additional species, including canine and bovine.

Suitable for use as control antibody for PARN siRNA (h): sc-61297, PARN siRNA (m): sc-61298, PARN shRNA Plasmid (h): sc-61297-SH, PARN shRNA Plasmid (m): sc-61298-SH, PARN shRNA (h) Lentiviral Particles: sc-61297-V and PARN shRNA (m) Lentiviral Particles: sc-61298-V.

PARN (N-20) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of PARN: 74 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.