

# BCAR3 (N-15): sc-47811

## BACKGROUND

BCAR3 (breast cancer anti-estrogen resistance 3, also designated AND-34 in murine; a novel SH2-containing protein 2 or NSP2, and SH2D3B) overexpression is a characteristic of antiestrogen resistance in human ZR-75-1 breast cancer cells. The deduced 825-amino acid BCAR3 protein contains a Src homology 2 (SH2) domain and shares homology to yeast CDC48. A 3.4-kb BCAR3 transcript is present in heart, placenta, skeletal muscle, spleen, prostate, testis, ovary, small intestine, colon, fetal kidney, and tamoxifen-resistant cancer cell lines. BCAR3 acts as a regulator of R-Ras to mediate the level of Insulin receptor substrate 1 (IRS-1) in MCF-7 and ZR-75-1 breast cancer cell lines. BCAR3 also interacts with p130Cas to enhance Src activation and cell migration.

## REFERENCES

1. Van Agthoven, T., et al. 1998. Identification of BCAR3 by a random search for genes involved in anti-estrogen resistance of human breast cancer cells. *EMBO J.* 17: 2799-2808.
2. Gotoh, T., et al. 2000. p130Cas regulates nucleotide exchange factor. *J. Biol. Chem.* 275: 30118-30123.

## CHROMOSOMAL LOCATION

Genetic locus: BCAR3 (human) mapping to 1p22.1; Bcar3 (mouse) mapping to 3 G1.

## SOURCE

BCAR3 (N-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of BCAR3 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-47811 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

BCAR3 (N-15) is recommended for detection of BCAR3 of mouse and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000). BCAR3 (N-15) is also recommended for detection of BCAR3 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for BCAR3 siRNA (h): sc-60265, BCAR3 siRNA (m): sc-60266, BCAR3 shRNA Plasmid (h): sc-60265-SH, BCAR3 shRNA Plasmid (m): sc-60266-SH, BCAR3 shRNA (h) Lentiviral Particles: sc-60265-V and BCAR3 shRNA (m) Lentiviral Particles: sc-60266-V.


Molecular Weight of BCAR3: 93 kDa.

Positive Controls: SK-BR-3 cell lysate: sc-2218, BT-20 cell lysate: sc-2223 or BCAR3 (h): 293T Lysate: sc-115467.


## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## DATA



BCAR3 (N-15): sc-47811. Western blot analysis of BCAR3 expression in non-transfected: sc-117752 (A) and human BCAR3 transfected: sc-115467 (B) 293T whole cell lysates.



BCAR3 (N-15): sc-47811. Immunofluorescence staining of methanol-fixed HeLa cells showing membrane and cytoplasmic localization.

## SELECT PRODUCT CITATIONS

1. Pennekamp, P., et al. 2011. Bcar3 is expressed in sertoli cells and germ cells of the developing testis in mice. *Sex. Dev.* 5: 197-204.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.

**MONOS**  
Satisfaction  
Guaranteed

Try **BCAR3 (3G4): sc-293346**, our highly recommended monoclonal alternative to BCAR3 (N-15).