

## FTCD (N-15): sc-47884

### BACKGROUND

58K protein antibodies are excellent for use as markers for the Golgi complex. The 58K protein has been identified as being FTCD, a bifunctional enzyme that channels 1-carbon units from formiminoglutamate, a metabolite of the histidine degradation pathway, to the folate pool. Defects in FTCD are the cause of glutamate formiminotransferase deficiency [also known as formiminoglutamicaciduria (FIGLU-uria)], an autosomal recessive disorder. Features of a severe phenotype include elevated levels of formiminoglutamate (FIGLU) in the urine in response to histidine administration, megaloblastic anemia and mental retardation. Features of a mild phenotype include high urinary excretion of FIGLU in the absence of histidine administration, mild developmental delay and no hematological abnormalities.

### REFERENCES

- Hennig, D., Scales, S.J., Moreau, A., Murley, L.L., De Mey, J. and Kreis, T.E. 1998. A formiminotransferase cyclodeaminase isoform is localized to the Golgi compl and can mediate interaction of *trans*-Golgi network-derived vesicles with microtubules. *J. Biol. Chem.* 273: 19602-19611.
- Bashour, A.M. and Bloom, G.S. 1998. 58K, a microtubule-binding Golgi protein, is a formiminotransferase cyclodeaminase. *J. Biol. Chem.* 273: 19612-19617.
- Gao, Y.S., Alvarez, C., Nelson, D.S. and Sztul, E. 1999. Molecular cloning, chara cyclodeaminase, a Golgi-associated 58 kDa protein. *J. Biol. Chem.* 273: 33825-33834.
- Gao, Y.S., Vrieling, A., MacKenzie, R. and Sztul, E. 2002. A novel type of regulation of the vimentin intermediate filament cytoskeleton by a Golgi protein. *Eur. J. Cell Biol.* 81: 391-401.
- Renous, R., Lapierre, P., Djilali-Saiah, I., Vitozzi, S. and Alvarez, F. 2003. In type 2 autoimmune hepatitis. *Exp. Cell Res.* 292: 332-341.
- Mao, Y., Vyas, N.K., Vyas, M.N., Chen, D.H., Ludtke, S.J., Chiu, W. and Quiocho, F.A. 2004. Structure of the bifunctional and Golgi-associated formiminotransferase octamer. *EMBO J.* 23: 2963-2971.
- Hagiwara, H., Tajika, Y., Matsuzaki, T., Suzuki, T., Aoki, T. and Takata, K. 2006. Localization of Golgi 58K protein (formiminotransferase cyclodeaminase) to the centrosome. *Histochem. Cell Biol.* 126: 251-259.

### CHROMOSOMAL LOCATION

Genetic locus: FTCD (human) mapping to 21q22.3; Ftcd (mouse) mapping to 10 C1.

### SOURCE

FTCD (N-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of FTCD of human origin.

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-47884 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

### APPLICATIONS

FTCD (N-15) is recommended for detection of FTCD isoforms A, D and C of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

FTCD (N-15) is also recommended for detection of FTCD isoforms A, D and C in additional species, including equine, bovine and porcine.

Suitable for use as control antibody for FTCD siRNA (h): sc-60662, FTCD siRNA (m): sc-60663, FTCD shRNA Plasmid (h): sc-60662-SH, FTCD shRNA Plasmid (m): sc-60663-SH, FTCD shRNA (h) Lentiviral Particles: sc-60662-V and FTCD shRNA (m) Lentiviral Particles: sc-60663-V.

Molecular Weight of FTCD: 58 kDa.

Positive Controls: rat liver extract: sc-2395.

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.


 MONOS  
Satisfation  
Guaranteed

Try **FTCD (G-3): sc-271788** or **FTCD (58K-9): sc-53128**, our highly recommended monoclonal alternatives to FTCD (N-15).