

# NDEL1 (C-20): sc-48063

## BACKGROUND

NudE-like protein (Ndel1) is expressed in the testis, brain, heart, hypothalamus, liver, lung, spleen, and stomach, specifically in the interphase centrosome and mitotic spindle. It positively regulates minus-end directed dynein. Evidence suggests that Ndel1 interacts with LIS1 to sustain the function of dynein, thereby impacting microtubule organization, nuclear translocation, and neuronal positioning. Ndel1 is phosphorylated during mitosis and seems to tether dynactin and dynein to the mother centriole for microtubule anchoring. Loss of function of Ndel1 in the developing neocortex impairs neuronal positioning and uncouples the centrosome and nucleus. Ndel1 may also impair mitochondrial transport or function, initiating a cascade of events culminating in psychiatric illness such as lissencephaly and schizophrenia.

## REFERENCES

- Niethammer, M., Smith, D.S., Ayala, R., Peng, J., Ko, J., Lee, M.S., Morabito, M. and Tsai, L.H. 2001. NUDEL is a novel Cdk5 substrate that associates with LIS1 and cytoplasmic dynein. *Neuron* 28: 697-711.
- Yan, X., Li, F., Liang, Y., Shen, Y., Zhao, X., Huang, Q. and Zhu, X. 2003. Human Nudel and NudE as regulators of cytoplasmic dynein in poleward protein transport along the mitotic spindle. *Mol. Cell. Biol.* 23: 1239-1250.
- Shu, T., Ayala, R., Nguyen, M.D., Xie, Z., Gleeson, J.G. and Tsai, L.H. 2004. NDEL1 operates in a common pathway with LIS1 and cytoplasmic dynein to regulate cortical neuronal positioning. *Neuron* 44: 263-277.
- Brandon, N.J., Schurov, I., Camargo, L.M., Handford, E.J., Duran-Jimeniz, B., Hunt, P., Millar, J.K., Porteous, D.J., Shearman, M.S. and Whiting, P.J. 2005. Subcellular targeting of DISC1 is dependent on a domain independent from the Nudel binding site. *Mol. Cell. Neurosci.* 28: 613-624.
- Li, J., Lee, W.L. and Cooper, J.A. 2005. NudEL targets dynein to microtubule ends through LIS1. *Nat. Cell Biol.* 7: 686-690.
- Toyo-Oka, K., Sasaki, S., Yano, Y., Mori, D., Kobayashi, T., Toyoshima, Y.Y., Tokuoka, S.M., Ishii, S., Shimizu, T., Muramatsu, M., Hiraiwa, N., Yoshiki, A., Wynshaw-Boris, A. and Hirotsune, S. 2005. Recruitment of katanin p60 by phosphorylated NDEL1, an LIS1 interacting protein, is essential for mitotic cell division and neuronal migration. *Hum. Mol. Genet.* 14: 3113-3128.

## CHROMOSOMAL LOCATION

Genetic locus: NDEL1 (human) mapping to 17p13.1; Ndel1 (mouse) mapping to 11 B3.

## SOURCE

NDEL1 (C-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of NDEL1 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-48063 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

NDEL1 (C-20) is recommended for detection of NDEL1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

NDEL1 (C-20) is also recommended for detection of NDEL1 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for NDEL1 siRNA (h): sc-61162, NDEL1 siRNA (m): sc-61163, NDEL1 shRNA Plasmid (h): sc-61162-SH, NDEL1 shRNA Plasmid (m): sc-61163-SH, NDEL1 shRNA (h) Lentiviral Particles: sc-61162-V and NDEL1 shRNA (m) Lentiviral Particles: sc-61163-V.

Molecular Weight of NDEL1: 40 kDa.

Positive Controls: HeLa + nocodazole cell lysate: sc-2274, MCF7 whole cell lysate: sc-2206 or SK-N-SH + nocodazole cell lysate: sc-2289.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.



Try **NDEL1 (D-5): sc-365094**, our highly recommended monoclonal alternative to NDEL1 (C-20).