

SBF2 (K-16): sc-48302

BACKGROUND

SET binding factor 2 (SBF2) is a pseudophosphatase that co-localizes with the myotubularin-related protein-2 (Mtmr2) forming a tetrameric complex in the cytoplasm of various tissues. SBF2 most likely plays a role in cellular communication or the signaling network that is necessary for myelin production, as well as in the development of the meshlike drainage canals surrounding the iris. Mutations the SBF2 gene have been identified as cause of charcot-marie-tooth disease type 4B2 (CMT4B2), an autosomal recessive demyelinating disease. Patients with CMT4B2 due to nonsense or truncating mutations in the SBF2 gene tend to develop early-onset open-angle glaucoma because of the complete absence of SBF2.

REFERENCES

1. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607697. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
2. Senderek, J., Bergmann, C., Weber, S., Ketelsen, U.P., Schorle, H., Rudnik-Schöneborn, S., Büttner, R., Buchheim, E. and Zerres, K. 2003. Mutation of the SBF2 gene, encoding a novel member of the myotubularin family, in Charcot-Marie-Tooth neuropathy type 4B2/11p15. *Hum. Mol. Genet.* 12: 349-356.
3. Hirano, R., Takashima, H., Umehara, F., Arimura, H., Michizono, K., Okamoto, Y., Nakagawa, M., Boerkoel, C.F., Lupski, J.R., Osame, M. and Arimura, K. 2004. SET binding factor 2 (SBF2) mutation causes CMT4B with juvenile onset glaucoma. *Neurology* 63: 577-580.
4. Conforti, F.L., Muglia, M., Mazzei, R., Patitucci, A., Valentino, P., Bono, F., Magariello, A., Sprovieri, T., Bergmann, C., Gabriele, A.L., Peluso, G., Nisticò, R., Senderek, J. and Quattrone, A. 2004. A new SBF2 mutation in a family with recessive demyelinating Charcot-Marie-Tooth (CMT4B2). *Neurology* 63: 1327-1328.
5. Srivastava, S., Li, Z., Lin, L., Liu, G., Ko, K., Coetzee, W.A. and Skolnik, E.Y. 2005. The phosphatidylinositol 3-phosphate phosphatase myotubularin-related protein 6 (MTMR6) is a negative regulator of the Ca²⁺-activated K⁺ channel KCa3.1. *Mol. Cell. Biol.* 25: 3630-3638.
6. Robinson, F.L. and Dixon, J.E. 2005. The phosphoinositide-3-phosphatase MTMR2 associates with MTMR13, a membrane-associated pseudophosphatase also mutated in type 4B Charcot-Marie-Tooth disease. *J. Biol. Chem.* 280: 31699-31707.
7. Berger, P., Berger, I., Schaffitzel, C., Tersar, K., Volkmer, B. and Suter, U. 2006. Multi-level regulation of myotubularin-related protein-2 phosphatase activity by myotubularin-related protein-13/SET-binding factor 2. *Hum. Mol. Genet.* 15: 569-579.

CHROMOSOMAL LOCATION

Genetic locus: SBF2 (human) mapping to 11p15.4; Sbf2 (mouse) mapping to 7 F2.

SOURCE

SBF2 (K-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of SBF2 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-48302 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

SBF2 (K-16) is recommended for detection of SBF2 of mouse and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

SBF2 (K-16) is also recommended for detection of SBF2 in additional species, including equine and porcine.

Suitable for use as control antibody for SBF2 siRNA (h): sc-61495, SBF2 siRNA (m): sc-61496, SBF2 shRNA Plasmid (h): sc-61495-SH, SBF2 shRNA Plasmid (m): sc-61496-SH, SBF2 shRNA (h) Lentiviral Particles: sc-61495-V and SBF2 shRNA (m) Lentiviral Particles: sc-61496-V.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.