

GFAT1 (C-14): sc-48512

BACKGROUND

Glutamine:fructose-6-phosphate amidotransferase (GFAT1) is the first and rate-limiting enzyme for the entry of glucose into the hexosamine biosynthesis pathway (HBP) in mammals. GFAT1, a member of the N-terminal nucleophile class of amidotransferases, converts fructose-6-phosphate into N-acetylglucosamine-6-phosphate. Hyperglycemia-induced Insulin resistance, a condition in which exposure to high concentrations of glucose and Insulin results in Insulin resistance, may result from increased glucose metabolism through the HBP. Hyperglycemia-induced Insulin resistance is a characteristic feature of type 2 diabetes. Consequently, GFAT1 is a potential therapeutic target in the treatment of type 2 diabetes.

REFERENCES

1. Nerlich, A.G., Sauer, U., Kolm-Litty, V., Wagner, E., Koch, M. and Schleicher, E.D. 1998. Expression of glutamine:fructose-6-phosphate amidotransferase in human evidence for high variability and distinct regulation in diabetes. *Diabetes* 47: 170-178.
2. DeHaven, J.E., Robinson, K.A., Nelson, B.A. and Buse, M.G. 2001. A novel variant of mRNA is selectively expressed in striated muscle. *Diabetes* 50: 2419-2424.
3. Niimi, M., Ogawara, T., Yamashita, T., Yamamoto, Y., Ueyama, A., Kambe, T., Okamoto, T., Ban, T., Tamanoi, H., Ozaki, K., Fujiwara, T., Fukui, H., Takahashi, E.I., Kyushiki, H. and Tanigami, A. 2001. Identification of GFAT1-L, a novel splice variant of human glutamine: fructose-6-phosphate amidotransferase (GFAT1) that is expressed abundantly in skeletal muscle. *J. Hum. Genet.* 46: 566-571.
4. Kaneto, H., Xu, G., Song, K.H., Suzuma, K., Bonner-Weir, S., Sharma, A. and Weir, G.C. 2001. Activation of the hexosamine pathway leads to deterioration of pancreatic β -cell function through the induction of oxidative stress. *J. Biol. Chem.* 276: 31099-31104.
5. Broschat, K.O., Gorka, C., Kasten, T.P., Gulve, EA. and Kilpatrick, B. 2002. A radiometric assay for glutamine:fructose-6-phosphate amidotransferase. *Anal. Biochem.* 305: 10-15.
6. Broschat, K.O., Gorka, C., Page, J.D., Martin-Berger, C.L., Davies, M.S., Huang Hc, H.C., Gulve, E.A., Salsgiver, W.J. and Kasten, T.P. 2002. Kinetic characterization of human glutamine-fructose-6-phosphate amidotransferase I: potent feedback inhibition by glucosamine 6-phosphate. *J. Biol. Chem.* 277: 14764-14770.
7. Chou, K.C. 2004. Molecular therapeutic target for type-2 diabetes. *J. Proteome Res.* 3: 1284-1288.
8. Elbein, S.C. Zheng, H., Jia, Y., Chu, W., Cooper, J.J., Hale, T. and Zhang, Z. 2004. Molecular screening of the human glutamine-fructose-6-phosphate amidotransferase 1 (GFPT1) gene and association studies with diabetes and diabetic nephropathy. *Mol. Genet. Metab.* 82: 321-328.

CHROMOSOMAL LOCATION

Genetic locus: GFPT1 (human) mapping to 2p13.3; Gfpt1 (mouse) mapping to 6 D1.

SOURCE

GFAT1 (C-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of GFAT1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-48512 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

GFAT1 (C-14) is recommended for detection of GFAT1 isoforms 1 and 2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

GFAT1 (C-14) is also recommended for detection of GFAT1 isoforms 1 and 2 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for GFAT1 siRNA (h): sc-60681, GFAT1 siRNA (m): sc-60682, GFAT1 shRNA Plasmid (h): sc-60681-SH, GFAT1 shRNA Plasmid (m): sc-60682-SH, GFAT1 shRNA (h) Lentiviral Particles: sc-60681-V and GFAT1 shRNA (m) Lentiviral Particles: sc-60682-V.

Molecular Weight of GFAT1: 77 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **GFAT1 (D-9): sc-377479**, our highly recommended monoclonal alternative to GFAT1 (C-14).