

Lsk/Ctk (H-87): sc-48827

BACKGROUND

All members of the Src gene family of tyrosine kinases are characterized by a carboxy terminal domain tyrosine, Y527 in the case of Src p60, which is highly phosphorylated in the inactive form of the enzyme, while phosphorylated to a much lesser extent when the enzyme is active. For instance, a mutant of c-Src, in which Y527 is replaced by phenylalanine, is transforming and displays 5 to 10-fold elevated kinase activity compared to its normal counterpart. Csk has been identified as a Src related tyrosine kinase having both SH2 and SH3 domains and a catalytic domain but lacking sequences amino terminal to the SH3 domain as well as the carboxy-terminal regulatory sequences. Csk phosphorylates Src on Y527 and also down regulates Lyn, Fyn and Lck by tyrosine phosphorylation of carboxy-terminal regulatory sites. A Csk-like protein-tyrosine kinase of mouse origin (Ctk), also designated Ntk, and its human homolog, Lsk, have also been described.

REFERENCES

- Okada, M., et al. 1989. A protein tyrosine kinase involved in regulation of pp60c-src function. *J. Biol. Chem.* 264: 20886-20893.
- Nada, S., et al. 1991. Cloning of a complementary DNA for a protein-tyrosine kinase that specifically phosphorylates a negative regulatory site of p60c-src. *Nature* 351: 69-72.
- Cooper, J.A., et al. 1993. The when and how of Src regulation. *Cell* 73: 1051-1054.
- Imamoto, A., et al. 1993. Disruption of the csk gene, encoding a negative regulator of Src family tyrosine kinases, leads to neural tube defects and embryonic lethality in mice. *Cell* 73: 1117-1124.
- Nada, S., et al. 1993. Constitutive activation of Src family kinases in mouse embryos that lack csk. *Cell* 73: 1125-1135.
- Superti-Furga, G., et al. 1993. Csk inhibition of c-Src activity requires both the SH2 and SH3 domains of Src. *EMBO J.* 12: 2625-2634.
- Klages, S., et al. 1994. Ctk: a protein-tyrosine kinase related to Csk that defines an enzyme family. *Proc. Natl. Acad. Sci. USA* 91: 2597-2601.
- Chow, L.M.L., et al. 1994. Ntk: a csk-related protein-tyrosine kinase expressed in brain and T lymphocytes. *Proc. Natl. Acad. Sci. USA* 91: 4975-4979.
- McVicar, D.W., et al. 1994. Molecular cloning of Lsk, a carboxyl-terminal Src kinase (csk) related gene, expressed in leukocytes. *Oncogene* 9: 2037-2044.

CHROMOSOMAL LOCATION

Genetic locus: MATK (human) mapping to 19p13.3; Matk (mouse) mapping to 10 C1.

SOURCE

Lsk/Ctk (H-87) is a rabbit polyclonal antibody raised against amino acids 171-257 mapping within an internal region of Lsk of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

Lsk/Ctk (H-87) is recommended for detection of Lsk of human origin and Ctk of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Lsk/Ctk (M-83) is also recommended for detection of Lsk in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for Lsk siRNA (h): sc-38973, Ctk siRNA (m): sc-38972, Lsk shRNA Plasmid (h): sc-38973-SH, Ctk shRNA Plasmid (m): sc-38972-SH, Lsk shRNA (h) Lentiviral Particles: sc-38973-V and Ctk shRNA (m) Lentiviral Particles: sc-38972-V.

Molecular Weight of Lsk/Ctk: 57 kDa.

Positive Controls: AML-193 whole cell lysate: sc-364182.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
 Satisfation
 Guaranteed

Try **Lsk/Ctk (H-1): sc-271174** or **Lsk/Ctk (32): sc-136309**, our highly recommended monoclonal alternatives to Lsk/Ctk (H-87).