

SNX17 (P-14): sc-49128

BACKGROUND

Sorting nexin (SNX) proteins are members of a large family of hydrophilic proteins that interact with a variety of receptor types, are involved in intracellular trafficking and contain a characteristic phox homology (PX) domain. SNX17, which demonstrates ubiquitous expression, contains a PX domain that shares 28% sequence identity with the PX domain of SNX1, as well as a B41 (FERM) domain. The SNX17 gene maps to chromosome 2 and is part of the cellular sorting machinery that regulates cell surface levels of LRP (lipoprotein receptor-related protein) by promoting its recycling. While the PX domain of SNX17 interacts with phosphatidylinositol-3-phosphate for membrane association, the FERM domain and the carboxyl-terminal region aid in LRP binding. Research indicates that SNX17 is localized to the limiting membrane and recycling tubules of early endosomes.

REFERENCES

- Nomura, N., Nagase, T., Miyajima, N., Sazuka, T., Tanaka, A., Sato, S., Seki, N., Kawarabayasi, Y., Ishikawa, K. and Tabata, S. 1994. Prediction of the coding sequences of unidentified human genes. II. The coding sequences of 40 new genes (KIAA0041-KIAA0080) deduced by analysis of cDNA clones from human cell line KG-1. *DNA Res.* 1: 223-229.
- Florian, V., Schlüter, T. and Bohnensack, R. 2001. A new member of the sorting nexin family interacts with the C-terminus of P-Selectin. *Biochem. Biophys. Res. Commun.* 281: 1045-1050.
- Stockinger, W., Sailler, B., Strasser, V., Recheis, B., Fasching, D., Kahr, L., Schneider, W.J. and Nimpf, J. 2002. The PX domain protein SNX17 interacts with members of the LDL receptor family and modulates endocytosis of the LDL receptor. *EMBO J.* 21: 4259-4267.
- Burden, J.J., Sun, X.M., García, A.B. and Soutar, A.K. 2004. Sorting motifs in the intracellular domain of the low density lipoprotein receptor interact with a novel domain of sorting nexin 17. *J. Biol. Chem.* 279: 16237-16245.
- Williams, R., Schlüter, T., Roberts, M.S., Knauth, P., Bohnensack, R. and Cutler, D.F. 2004. Sorting nexin 17 accelerates internalization yet retards degradation of P-Selectin. *Mol. Biol. Cell* 15: 3095-3105.
- Knauth, P., Schlüter, T., Czubayko, M., Kirsch, C., Florian, V., Schreckenberger, S., Hahn, H. and Bohnensack, R. 2005. Functions of sorting nexin 17 do-mains and recognition motif for P-Selectin trafficking. *J. Mol. Biol.* 347: 813-825.
- van Kerkhof, P., Lee, J., McCormick, L., Tetrault, E., Lu, W., Schoenfish, M., Oorschot, V., Strous, G.J., Klumperman, J. and Bu, G. 2005. Sorting nexin 17 facilitates LRP recycling in the early endosome. *EMBO J.* 24: 2851-2861.

CHROMOSOMAL LOCATION

Genetic locus: SNX17 (human) mapping to 2p23-p22; Snx17 (mouse) mapping to 5 B1.

SOURCE

SNX17 (P-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of SNX17 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-49128 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

SNX17 (P-14) is recommended for detection of SNX17 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


Suitable for use as control antibody for SNX17 siRNA (h): sc-61587, SNX17 siRNA (m): sc-61588, SNX17 shRNA Plasmid (h): sc-61587-SH, SNX17 shRNA Plasmid (m): sc-61588-SH, SNX17 shRNA (h) Lentiviral Particles: sc-61587-V and SNX17 shRNA (m) Lentiviral Particles: sc-61588-V.

Molecular Weight of SNX17: 52.8 kDa.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible mouse anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



SNX17 (P-14): sc-49128. Western blot analysis of SNX17 expression in non-transfected: sc-117752 (A) and mouse SNX17 transfected: sc-123694 (B) 293T whole cell lysates.



SNX17 (P-14): sc-49128. Western blot analysis of SNX17 expression in non-transfected: sc-117752 (A) and human SNX17 transfected: sc-116283 (B) 293T whole cell lysates.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.