

UROS (C-16): sc-49568

BACKGROUND

Uroporphyrinogen-III synthase (UROS), also designated hydroxymethylbilane hydrolyase, is a 265-amino acid ubiquitously expressed protein that is predominantly found in skeletal muscle, heart muscle, the caudate nucleus, and the amygdala. UROS is the fourth enzyme in the human heme biosynthetic pathway and is responsible for the conversion of hydroxymethylbilane to the cyclic tetrapyrrole, uroporphyrinogen III. UROS is a monomer involved in porphyrin biosynthesis by the Shemin pathway. Defects in the UROS gene can cause molecular lesions that lead to the autosomal recessive Gunther disease, also known as congenital erythropoietic porphyria (CEP). CEP is characterized by accumulation of uroporphyrin I in the bone marrow, peripheral blood, and other organs, often causing death early in adult life.

REFERENCES

1. Tsai, S.F., Bishop, D.F. and Desnick, R.J. 1987. Purification and properties of uroporphyrinogen III synthase from human erythrocytes. *J. Biol. Chem.* 262: 1268-1273.
2. Tsai, S.F., Bishop, D.F. and Desnick, R.J. 1988. Human uroporphyrinogen III synthase: molecular cloning, nucleotide sequence, and expression of a full-length cDNA. *Proc. Natl. Acad. Sci. USA* 85: 7049-7053.
3. Astrin, K.H., Warner, C.A., Yoo, H.W., Goodfellow, P.J., Tsai, S.F. and Desnick, R.J. 1991. Regional assignment of the to chromosome 10q25.2-q26.3. *Hum. Genet.* 87: 18-22.
4. Xu, W., Kozak, C.A. and Desnick, R.J. 1995. Uroporphyrinogen-III synthase: molecular cloning, nucleotide sequence, expression of a mouse full-length cDNA, and its localization on mouse chromosome 7. *Genomics* 26: 556-562.
5. Geronimi, F., Richard, E., Lamrissi-Garcia, I., Lalanne, M., Ged, C., Redonnet-Vernhet, I., Moreau-Gaudry, F. and de Verneuil, H. 2003. Lentivirus-mediated gene transfer of uroporphyrinogen III synthase fully corrects the porphyric phenotype in human cells. *J. Mol. Med.* 81: 310-320.
6. Torres, G., Quiñones, M., Hau, J., Perez, L. and García, O. 2003. HLA-DQA1 and polymarker in the UROS population from Peruvian Highland. *J. Forensic Sci.* 48: 1426-1427.
7. Omata, Y., Sakamoto, H., Higashimoto, Y., Hayashi, S. and Noguchi, M. 2004. Purification and characterization of human uroporphyrinogen III synthase expressed in *Escherichia coli*. *J. Biochem.* 136: 211-220.
8. de Villiers, J.N., Kotze, M.J., van Heerden, C.J., Sadie, A., Gardner, H.F., Liebenberg, J., van Zyl, R., du Plessis, L., Kimberg, M., Frank, J. and Warnich, L. 2005. Overrepresentation of the founder PPOX gene mutation R59W in a South African patient with severe clinical manifestation of porphyria. *Exp. Dermatol.* 14: 50-55.

CHROMOSOMAL LOCATION

Genetic locus: UROS (human) mapping to 10q26.13; Uros (mouse) mapping to 7 F3.

SOURCE

UROS (C-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of UROS of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-49568 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

UROS (C-16) is recommended for detection of Uroporphyrinogen-III synthase (UROS) of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

UROS (C-16) is also recommended for detection of Uroporphyrinogen-III synthase (UROS) in additional species, including equine, canine and porcine.

Suitable for use as control antibody for UROS siRNA (h): sc-61760, UROS siRNA (m): sc-61761, UROS shRNA Plasmid (h): sc-61760-SH, UROS shRNA Plasmid (m): sc-61761-SH, UROS shRNA (h) Lentiviral Particles: sc-61760-V and UROS shRNA (m) Lentiviral Particles: sc-61761-V.

Molecular Weight of UROS: 30 kDa.

Positive Controls: HL-60 whole cell lysate: sc-2209, mouse liver extract: sc-2256 or mouse thyroid extract: sc-2407.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.