

IRSp53 (C-16): sc-50008

BACKGROUND

The scaffolding protein Insulin receptor tyrosine kinase substrate p53 (IRSp53), an ubiquitous regulator of the Actin cytoskeleton, mediates filopodia formation under the control of Rho-family GTPases. It is expressed in the cytoplasm and links small membrane-bound G proteins to cytoplasmic effector proteins. IRSp53 comprises a central SH3 domain, which binds to proline-rich regions of a wide range of Actin regulators, and a conserved N-terminal IRSp53/MIM homology domain (IMD) that harbors F-Actin-bundling activity. IRSp53 interacts with atrophin-1, the product of the dentatorubral-pallidoluysian atrophy (DRPLA) gene, which is associated with an autosomal dominant neurodegenerative disease. The IRSp53 protein also interacts with ENAH, BAI1, Eps8, Shank 1, Shank 2, Shank 3, WAVE1, WAVE2, Tiam1 and Dia 1.

REFERENCES

- Okamura-Oho, Y., et al. 1999. Dentatorubral-pallidoluysian atrophy protein interacts through a proline-rich region near polyglutamine with the SH3 domain of an Insulin receptor tyrosine kinase substrate. *Hum. Mol. Genet.* 8: 947-957.
- Soltau, M., et al. 2002. The Insulin receptor substrate IRSp53 links post-synaptic Shank 1 to the small G protein Cdc42. *Mol. Cell. Neurosci.* 21: 575-583.
- Miyahara, A., et al. 2003. Genomic structure and alternative splicing of the Insulin receptor tyrosine kinase substrate of 53 kDa protein. *J. Hum. Genet.* 48: 410-414.
- Funato, Y., et al. 2004. IRSp53/Eps8 complex is important for positive regulation of Rac and cancer cell motility/invasiveness. *Cancer Res.* 64: 5237-5244.
- Choi, J. et al. 2005. Regulation of dendritic spine morphogenesis by Insulin receptor substrate 53, a downstream effector of Rac1 and Cdc42 small GTPases. *J. Neurosci.* 25: 869-789.
- Connolly, B.A. et al. 2005. Tiam1-IRSp53 complex formation directs specificity of Rac-mediated Actin cytoskeleton regulation. *Mol. Cell. Biol.* 25: 4602-4614.
- Hori, K. et al. 2005. NMDA receptor-dependent synaptic translocation of Insulin receptor substrate p53 via protein kinase C signaling. *J. Neurosci.* 25: 2670-2681.
- Mackie, S. et al. 2005. Novel brain 14-3-3 interacting proteins involved in neurodegenerative disease. *FEBS J.* 272: 4202-4210.

CHROMOSOMAL LOCATION

Genetic locus: BAIAP2 (human) mapping to 17q25.3; Baiap2 (mouse) mapping to 11 E2.

SOURCE

IRSp53 (C-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of IRSp53 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-50008 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

IRSp53 (C-16) is recommended for detection of IRSp53 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

IRSp53 (C-16) is also recommended for detection of IRSp53 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for IRSp53 siRNA (h): sc-60863, IRSp53 siRNA (m): sc-60864, IRSp53 shRNA Plasmid (h): sc-60863-SH, IRSp53 shRNA Plasmid (m): sc-60864-SH, IRSp53 shRNA (h) Lentiviral Particles: sc-60863-V and IRSp53 shRNA (m) Lentiviral Particles: sc-60864-V.

Molecular Weight of IRSp53: 53 kDa.

Positive Controls: mouse brain extract: sc-2253, K-562 whole cell lysate: sc-2203 or mouse testis extract: sc-2405.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

SELECT PRODUCT CITATIONS

- Romero, A.M., et al. 2010. Chronic ethanol exposure alters the levels, assembly, and cellular organization of the actin cytoskeleton and microtubules in hippocampal neurons in primary culture. *Toxicol. Sci.* 18: 602-612.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.



Try **IRSp53 (46): sc-136470**, our highly recommended monoclonal alternative to IRSp53 (C-16).