

Lipin-1 (H-19): sc-50050

BACKGROUND

The Lipin family of nuclear proteins contains three members: Lipin-1, Lipin-2 and Lipin-3, all of which contain a nuclear signal sequence, a highly conserved amino-terminal (NLIP) domain and a carboxy-terminal (CLIP) domain. Lipin-1 is crucial for normal adipose tissue development and metabolism. Lipin-1 selectively activates a subset of PGC-1 α target pathways, including fatty acid oxidation and mitochondrial oxidative phosphorylation, by inducing expression of the nuclear receptor PPAR α . Lipin-1 also inactivates the lipogenic program and suppresses circulating lipid levels. An abundance of Lipin-1 promotes fat accumulation and Insulin sensitivity, whereas a deficiency in Lipin-1 may deter normal adipose tissue development, resulting in Insulin resistance and lipodystrophy, a heterogeneous group of disorders characterized by loss of body fat, fatty liver, hypertriglyceridemia and Insulin resistance.

REFERENCES

- Peterfy, M., Phan, J., Xu, P. and Reue, K. 2001. Lipodystrophy in the fld mouse results from mutation of a new gene encoding a nuclear protein, Lipin. *Nat. Genet.* 27: 121-124.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 605518. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Reitman, M.L. 2005. The fat and thin of Lipin. *Cell Metab.* 1: 5-6.
- Phan, J. and Reue, K. 2005. Lipin, a lipodystrophy and obesity gene. *Cell Metab.* 1: 73-83.
- Phan, J., Peterfy, M. and Reue, K. 2005. Biphasic expression of Lipin suggests dual roles in adipocyte development. *Drug News Perspect.* 18: 5-11.
- Finck, B.N., Gropler, M.C., Chen, Z., Leone, T.C., Croce, M.A., Harris, T.E., Lawrence, J.C. and Kelly, D.P. 2006. Lipin-1 is an inducible amplifier of the hepatic PGC-1 α /PPAR α regulatory pathway. *Cell Metab.* 4: 199-210.
- Han, G.S., Wu, W.I. and Carman, G.M. 2006. The *Saccharomyces cerevisiae* Lipin homolog is a Mg²⁺-dependent phosphatidate phosphatase enzyme. *J. Biol. Chem.* 281: 9210-9218.
- Larosa, P.C., Miner, J., Xia, Y., Zhou, Y., Kachman, S. and Fromm, M.E. 2006. *Trans*-10, *cis*-12 conjugated linoleic acid causes inflammation and delipidation of white adipose tissue in mice: a microarray and histological analysis. *Physiol. Genomics* 27: 282-294.
- Suviolahti, E., Reue, K., Cantor, R.M., Phan, J., Gentile, M., Naukkarinen, J., Soro-Paavonen, A., Oksanen, L., Kaprio, J., Rissanen, A., Salomaa, V., Kontula, K., Taskinen, M.R., Pajukanta, P. and Peltonen, L. 2006. Cross-species analyses implicate Lipin-1 involvement in human glucose metabolism. *Hum. Mol. Genet.* 15: 377-386.

CHROMOSOMAL LOCATION

Genetic locus: LPIN1 (human) mapping to 2p25.1; Lpin1 (mouse) mapping to 12 A1.1.

SOURCE

Lipin-1 (H-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of Lipin-1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-50050 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Lipin-1 (H-19) is recommended for detection of Lipin-1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Lipin-1 (H-19) is also recommended for detection of Lipin-1 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for Lipin-1 siRNA (h): sc-60940, Lipin-1 siRNA (m): sc-60941, Lipin-1 shRNA Plasmid (h): sc-60940-SH, Lipin-1 shRNA Plasmid (m): sc-60941-SH, Lipin-1 shRNA (h) Lentiviral Particles: sc-60940-V and Lipin-1 shRNA (m) Lentiviral Particles: sc-60941-V.

Molecular Weight of Lipin-1: 102 kDa.

Positive Controls: Jurkat whole cell lysate: sc-2204.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **Lipin-1 (B-12): sc-376874**, our highly recommended monoclonal alternative to Lipin-1 (H-19).