SANTA CRUZ BIOTECHNOLOGY, INC.

AHSP (M-18): sc-50793



BACKGROUND

α-hemoglobin stabilizing protein (AHSP), also designated erythroid associated factor (ERAF), is an erythroid-specific protein that acts as a chaperone to prevent the aggregation of α-hemoglobin during normal erythroid cell development. It specifically protects free α-hemoglobin from precipitation in live cells and in solution. It forms a heterodimer with free α-hemoglobin, but not with β-hemoglobin or hemoglobin A (α2-β-2). AHSP localizes to the cytoplasm and is expressed in the blood and bone marrow. The AHSP protein is down-regulated in transmissible spongiform encephalopathies (TSEs). AHSP may regulate pathological states of α-hemoglobin excess such as β-thalassemia, a group of hereditary disorders involving the decreased production of normal adult hemoglobin (HbA) that are characterized by a deficiency in the synthesis of β-globin chains.

REFERENCES

- Zhang, Q.H., Ye, M., Wu, X.Y., Ren, S.X., Zhao, M., Zhao, C.J., Fu, G., Shen, Y., Fan, H.Y., Lu, G., Zhong, M., Xu, X.R., Han, Z.G., Zhang, J.W., Tao, J., Huang, Q.H., Zhou, J., Hu, G.X., Gu, J., Chen, S.J. and Chen, Z. 2000. Cloning and functional analysis of cDNAs with open reading frames for 300 previously undefined genes expressed in CD34⁺ hematopoietic stem/ progenitor cells. Genome Res. 10: 1546-1560.
- Miele, G., Manson, J. and Clinton, M. 2001. A novel erythroid-specific marker of transmissible spongiform encephalopathies. Nat. Med. 7: 361-364.
- Kihm, A.J., Kong, Y., Hong, W., Russell, J.E., Rouda, S., Adachi, K., Simon, M.C., Blobel, G.A. and Weiss, M.J. 2002. An abundant erythroid protein that stabilizes free α-haemoglobin. Nature 417: 758-763.
- Gell, D., Kong, Y., Eaton, S.A., Weiss, M.J. and Mackay, J.P. 2002. Biophysical characterization of the α-globin binding protein α-hemoglobin stabilizing protein. J. Biol. Chem. 277: 40602-40609.
- 5. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 605821. World Wide Web URL: http://www.ncbi.nlm.nih.gov/omim/
- dos Santos, C.O. and Costa, F.F. 2005. AHSP and β-thalassemia: a possible genetic modifier. Hematology 10: 157-161.

SOURCE

AHSP (M-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of AHSP of mouse origin.

PRODUCT

Each vial contains 200 μg lgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-50793 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

AHSP (M-18) is recommended for detection of AHSP of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for AHSP siRNA (m): sc-60138, AHSP shRNA Plasmid (m): sc-60138-SH and AHSP shRNA (m) Lentiviral Particles: sc-60138-V.

Molecular Weight of AHSP: 12 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS Satisfation Guaranteed

Try **AHSP (G-5): sc-515189**, our highly recommended monoclonal alternative to AHSP (M-18).