

CA V (R-18): sc-50818

BACKGROUND

Carbonic anhydrases (CAs) are members of a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. CAs are involved in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption and the formation of aqueous humor, cerebrospinal fluid, saliva and gastric juice. They show extensive diversity in distribution and in their subcellular localization. The human CA5A gene encodes for the CA V protein which is expressed in the mitochondrial matrix in liver tissues, but can also be detected in many other tissues. In ureagenesis, the CA V protein may be involved in providing bicarbonate to carbamoyl phosphate synthetase.

REFERENCES

1. Dodgson, S.J., Tashian, R.E., Gross, G. and Carter, N.D. 1991. The Carbonic Anhydrases: Cellular Physiology and Molecular Genetics. New York: Plenum.
2. Sly, W.S. and Hu, P.Y. 1995. Human carbonic anhydrases and carbonic anhydrase deficiencies. *Annu. Rev. Biochem.* 64: 375-401.
3. Nagao, Y., Batanian, J.R., Clemente, M.F. and Sly, W.S. 1996. Genomic organization of the human gene (CA5) and pseudogene for mitochondrial carbonic anhydrase V and their localization to chromosomes 16q and 16p. *Genomics* 28: 477-484.
4. Duda, D.M., Tu, C., Fisher, S.Z., An, H., Yoshioka, C., Govindasamy, L., Laipis, P.J., Agbandje-McKenna, M., Silverman, D.N. and McKenna, R. 2005. Human carbonic anhydrase III: structural and kinetic study of catalysis and proton transfer. *Biochemistry* 44: 10046-10053.
5. Ishii, Y., Akazawa, D., Aoki, Y., Yamada, H. and Oguri, K. 2005. Suppression of carbonic anhydrase III mRNA level by an aryl hydrocarbon receptor ligand in primary cultured hepatocytes of rat. *Biol. Pharm. Bull.* 28: 1087-1090.
6. Liu, Y.G., Xion, Y.Z. and Deng C.Y. 2005. Isolation, sequence analysis and expression profile of a novel swine gene differentially expressed in the Longissimus dorsi muscle tissues from LandraceXLarge white cross-combination. *Acta. Biochim. Biophys. Sin.* 37: 186-191.
7. Masquelet, A.C. 2005. The anatomy lesson of Dr. Tulp. *J. Hand Surg.* 30: 379-381.

CHROMOSOMAL LOCATION

Genetic locus: Car5a (rat) mapping to 19q12.

SOURCE

CA V (R-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of CA V of rat origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-50818 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

CA V (R-18) is recommended for detection of CA V (carbonic anhydrase) of rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


Molecular Weight of CA V: 35 kDa.

Positive Controls: rat liver extract: sc-2395.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



CA V (R-18): sc-50818. Western blot analysis of CA V expression in rat liver tissue extract.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

MONOS
Satisfaction
Guaranteed

Try **CA V (A-1): sc-398184**, our highly recommended monoclonal alternative to CA V (R-18).