

NAGAT (T-17): sc-50850



The Power to Question

BACKGROUND

Histo-blood group ABO system transferase (NAGAT) is a member of the glycosyltransferase 6 family of proteins and the basis of the ABO blood group system. The histo-blood group ABO involves three carbohydrate antigens: A, B and H; the NAGAT protein functions as the basis of this group. A, B and AB individuals express a glycosyltransferase activity converting the H antigen to the A antigen (by addition of UDP-GalNAc) or to the B antigen (by addition of UDP-Gal), while O individuals do not express this glycosyltransferase activity. The B phenotype of NAGAT differs from the A form by a few residue substitutions, whereas the O form is a result of a single base frame-shift deletion in the N-terminal extremity of the gene. The NAGAT protein localizes to the Golgi apparatus and its conserved DXD motif functions in cofactor binding.

REFERENCES

1. Nydegger, U., Mohacsi, P., Koestner, S., Kappeler, A., Schaffner, T. and Carrel, T. 2004. ABO histo-blood group system-incompatible allografting. *Int. Immunopharmacol.* 5: 147-153.
2. Hutson, A.M., Airaud, F., LePendou, J., Estes, M.K. and Atmar, R.L. 2005. Norwalk virus infection associates with secretor status genotyped from sera. *J. Med. Virol.* 77: 116-120.
3. Nydegger, U.E., Tevaearai, H., Berdat, P., Rieben, R., Carrel, T., Mohacsi, P. and Flegel, W.A. 2005. Histo-blood group antigens as allo- and autoantigens. *Annu. N.Y. Acad. Sci.* 1050: 40-51.
4. Chen, Y.L., Chen, J.C., Lin, T.M., Huang, T.J., Wang, S.T., Lee, M.F. and Wang, J.Y. 2005. ABO/secretor genetic complex is associated with the susceptibility of childhood asthma in Taiwan. *Clin. Exp. Allergy* 35: 926-932.
5. Fretz, R., Svoboda, P., Schorr, D., Tanner, M. and Baumgartner, A. 2005. Risk factors for infections with Norovirus gastrointestinal illness in Switzerland. *Eur. J. Clin. Microbiol. Infect. Dis.* 24: 256-261.
6. Kominato, Y., Hata, Y., Matsui, K. and Takizawa, H. 2005. Regulation of ABO gene expression. *Leg. Med.* 7: 263-265.
7. Tanaka, Y., Haga, H., Egawa, H., Okuno, T., Miyagawa-Hayashino, A., Tsuruyama, T., Kambe, M., Marusawa, H., Chiba, T. and Manabe, T. 2005. Intra-graft expression of recipient-type ABO blood group antigens: long-term follow-up and histological features after liver transplantation. *Liver Transpl.* 11: 547-554.
8. Huang, P., Farkas, T., Zhong, W., Tan, M., Thornton, S., Morrow, A.L. and Jiang, X. 2005. Norovirus and histo-blood group antigens: demonstration of a wide spectrum of strain specificities and classification of two major binding groups among multiple binding patterns. *J. Virol.* 79: 6714-6722.

CHROMOSOMAL LOCATION

Genetic locus: ABO (human) mapping to 9q34.2; Abo (mouse) mapping to 2 A3.

SOURCE

NAGAT (T-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of NAGAT of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-50850 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

NAGAT (T-17) is recommended for detection of NAGAT of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

NAGAT (T-17) is also recommended for detection of NAGAT in additional species, including canine.

Suitable for use as control antibody for NAGAT siRNA (h): sc-61138, NAGAT siRNA (m): sc-61139, NAGAT shRNA Plasmid (h): sc-61138-SH, NAGAT shRNA Plasmid (m): sc-61139-SH, NAGAT shRNA (h) Lentiviral Particles: sc-61138-V and NAGAT shRNA (m) Lentiviral Particles: sc-61139-V.

Molecular Weight of NAGAT: 41 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.