



# Cryptic (T-15): sc-51062

## BACKGROUND

Cryptic is a 224 amino acid member of the epidermal growth factor-Cripto-1/FRL1/cryptic family that consists of extracellular and cell-associated proteins which function as co-receptors for the transforming growth factor  $\beta$  subfamily of proteins. Cryptic is involved in the appropriate establishment of the left-right axis of the organism, and it may play a role in the development of mesoderm and/or neural patterning during gastrulation. Loss-of-function mutations in CFC1, the gene encoding for Cryptic, are associated with human left-right laterality defects. Defects in CFC1 also cause visceral heterotaxy (HTX2), an autosomal disease characterized by a variable group of congenital anomalies that include complex cardiac malformations such as conotruncal heart malformations that may lead to cardiac outflow tract defects, pulmonary atresia, double-outlet right ventricle, truncus arteriosus communis and aortic arch anomalies.

## REFERENCES

- Burdine, R.D. and Schier, A.F. 2000. Conserved and divergent mechanisms in left-right axis formation. *Genes Dev.* 14: 763-776.
- Shen, M.M. and Schier, A.F. 2000. The EGF-CFC gene family in vertebrate development. *Trends Genet.* 16: 303-309.
- Bamford, R.N., Roessler, E., Burdine, R.D., Saplakoglu, U., dela Cruz, J., Splitt, M., Goodship, J.A., Towbin, J., Bowers, P., Ferrero, G.B., Marino, B., Schier, A.F., Shen, M.M., Muenke, M. and Casey, B. 2000. Loss-of-function mutations in the EGF-CFC gene CFC1 are associated with human left-right laterality defects. *Nat. Genet.* 26: 365-369.
- Goldmuntz, E., Bamford, R., Karkera, J.D., dela Cruz, J., Roessler, E. and Muenke, M. 2002. CFC1 mutations in patients with transposition of the great arteries and double-outlet right ventricle. *Am. J. Hum. Genet.* 70: 776-780.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 605194. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Hu, X.F. and Xing, P.X. 2005. Cripto as a target for cancer immunotherapy. *Expert Opin. Ther. Targets* 9: 383-394.
- Hu, X.F. and Xing, P.X. 2005. Cripto monoclonal antibodies. *Drug News Perspect.* 18: 293-303.
- Strizzi, L., Bianco, C., Normanno, N. and Salomon, D. 2005. Cripto-1: a multifunctional modulator during embryogenesis and oncogenesis. *Oncogene* 24: 5731-5741.
- Baroni, M., Mazzola, G., Kaabache, T., Borgel, D., Gandrille, S., Viganò, D'Angelo, S., Marchetti, G., di lasio, M.G., Pinotti, M., D'Angelo, A. and Bernardi, F. 2006. Molecular bases of type II protein S deficiency: the I203-D204 deletion in the EGF4 domain alters GLA domain function. *J. Thromb. Haemost.* 4: 186-191.

## CHROMOSOMAL LOCATION

Genetic locus: CFC1 (human) mapping to 2q21.1; Cfc1 (mouse) mapping to 1 B.

## SOURCE

Cryptic (T-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Cryptic of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-51062 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

Cryptic (T-15) is recommended for detection of Cryptic of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

Suitable for use as control antibody for Cryptic siRNA (h): sc-60453.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.