

MFF (A-7): sc-515664

BACKGROUND

MFF (mitochondrial fission factor), also known as GL004, AD030 or AD033, is a 342 amino acid single-pass type IV membrane protein of the mitochondrial outer membrane that belongs to the tango11 family. Involved in mitochondrial and peroxisomal fission, MFF is abundantly expressed in stomach, heart, muscle, liver, brain and kidney. MFF exists as five alternatively spliced isoforms that are encoded by a gene that maps to human chromosome 2q36.3. As the second largest human chromosome, chromosome 2 consists of 237 million bases, encodes over 1,400 genes and makes up approximately 8% of the human genome. A number of genetic diseases are linked to genes on chromosome 2. Harlequin ichthyosis, a rare and morbid skin deformity, is associated with mutations in the ABCA12 gene. The lipid metabolic disorder sitosterolemia is associated with ABCG5 and ABCG8. An extremely rare recessive genetic disorder, Alström syndrome, is due to mutations in the ALMS1 gene.

REFERENCES

- Patel, S.B., et al. 1998. Mapping a gene involved in regulating dietary cholesterol absorption. The sitosterolemia locus is found at chromosome 2p21. *J. Clin. Invest.* 102: 1041-1044.
- Zumsteg, U., et al. 2000. Alstrom syndrome: confirmation of linkage to chromosome 2p12-13 and phenotypic heterogeneity in three affected sibs. *J. Med. Genet.* 37: E8.
- Shulenin, S., et al. 2001. An ATP-binding cassette gene (ABCG5) from the ABCG (white) gene subfamily maps to human chromosome 2p21 in the region of the Sitosterolemia locus. *Cytogenet. Cell Genet.* 92: 204-208.
- Hearn, T., et al. 2002. Mutation of ALMS1, a large gene with a tandem repeat encoding 47 amino acids, causes Alström syndrome. *Nat. Genet.* 31: 79-83.

CHROMOSOMAL LOCATION

Genetic locus: MFF (human) mapping to 2q36.3; Mff (mouse) mapping to 1 C5.

SOURCE

MFF (A-7) is a mouse monoclonal antibody specific for an epitope mapping between amino acids 119-138 within a cytoplasmic domain of MFF of human origin.

PRODUCT

Each vial contains 200 µg IgG_{2a} kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-515664 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% stabilizer protein).

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

MFF (A-7) is recommended for detection of MFF of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for MFF siRNA (h): sc-94736, MFF siRNA (m): sc-149404, MFF shRNA Plasmid (h): sc-94736-SH, MFF shRNA Plasmid (m): sc-149404-SH, MFF shRNA (h) Lentiviral Particles: sc-94736-V and MFF shRNA (m) Lentiviral Particles: sc-149404-V.

Molecular Weight (predicted) of MFF: 38 kDa.


Molecular Weight (observed) of MFF: 25-39 kDa.

Positive Controls: MFF (m): 293T Lysate: sc-121620, HeLa whole cell lysate: sc-2200 or NCI-H460 whole cell lysate: sc-364235.


RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850. 4) Immunohistochemistry: use m-IgGκ BP-HRP: sc-516102 with DAB, 50X: sc-24982 and Immunohistomount: sc-45086, or Organo/Limonene Mount: sc-45087.

DATA



MFF (C-11): sc-515664. Western blot analysis of MFF expression in non-transfected 293T: sc-117752 (A), mouse MFF transfected 293T: sc-121620 (B), HeLa (C), NCI-H460 (D) and IMR-32 (E) whole cell lysates.



MFF (A-7): sc-515664. Immunoperoxidase staining of formalin fixed, paraffin-embedded mouse testis (A) and rat testis (B) tissue showing cytoplasmic staining of cells in seminiferous ducts and Leydig cells. Blocked with 0.25X UltraCruz® Blocking Reagent: sc-516214. Detection reagents used: m-IgGκ BP-B: sc-516142 and ImmunoCruz® ABC Kit: sc-516216.

RESEARCH USE

For research use only, not for use in diagnostic procedures.