

# Peroxin 19 (2E4): sc-517578

## BACKGROUND

The covalent attachment of prenyl lipids, such as farnesyl or geranylgeranyl groups, by specific transferases is indispensable for the cellular sorting of many proteins. A farnesylated protein, Peroxin 19 (peroxisomal farnesylated protein, PxF or HK33), localizes to the outer surface of peroxisomes. Specifically, Peroxin 19 localizes to the cytoplasmic surface of peroxisomes in liver cells. Peroxin 19 is the human homolog of Pex19p, an oleic acid-inducible, farnesylated protein that is essential for peroxisome biogenesis in *Saccharomyces cerevisiae*. The carboxy-terminal part of Peroxin 19, including the CAAX homology box, is required for its biologic function. Moreover, Peroxin 19 is apparently involved in the initial stage of peroxisome membrane assembly, before the import of matrix protein. The gene which encodes Peroxin 19 is a housekeeping gene and maps to human chromosome 1q23.2. This is the causative gene for complementation group J peroxisome biogenesis disorder.

## REFERENCES

- James, G.L., Goldstein, J.L., Pathak, R.K., Anderson, R.G. and Brown, M.S. 1994. PxF, a prenylated protein of peroxisomes. *J. Biol. Chem.* 269: 14182-14190.
- Kammerer, S., Arnold, N., Gutensohn, W., Mewes, H.W., Kunau, W.H., Hofler, G., Roscher, A.A. and Braun, A. 1997. Genomic organization and molecular characterization of a gene encoding HsPXF, a human peroxisomal farnesylated protein. *Genomics* 45: 200-210.
- Gotte, K., Girzalsky, W., Linkert, M., Baumgart, E., Kammerer, S., Kunau, W.H. and Erdmann, R. 1998. Pex19p, a farnesylated protein essential for peroxisome biogenesis. *Mol. Cell. Biol.* 18: 616-628.
- Kinoshita, N., Ghaedi, K., Shimozawa, N., Wanders, R.J., Matsuzono, Y., Imanaka, T., Okumoto, K., Suzuki, Y., Kondo, N. and Fujiki, Y. 1998. Newly identified Chinese hamster ovary cell mutants are defective in biogenesis of peroxisomal membrane vesicles (Peroxisomal ghosts), representing a novel complementation group in mammals. *J. Biol. Chem.* 273: 24122-24130.
- Matsuzono, Y., Kinoshita, N., Tamura, S., Shimozawa, N., Hamasaki, M., Ghaedi, K., Wanders, R.J., Suzuki, Y., Kondo, N. and Fujiki, Y. 1999. Human PEX19: cDNA cloning by functional complementation, mutation analysis in a patient with Zellweger syndrome, and potential role in peroxisomal membrane assembly. *Proc. Natl. Acad. Sci. USA* 96: 2116-2121.

## CHROMOSOMAL LOCATION

Genetic locus: PEX19 (human) mapping to 1q23.2.

## SOURCE

Peroxin 19 (2E4) is a mouse monoclonal antibody raised against full length recombinant Peroxin 19 of human origin.

## PRODUCT

Each vial contains 100 µg IgG<sub>2a</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

Peroxin 19 (2E4) is recommended for detection of Peroxin 19 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Peroxin 19 siRNA (h): sc-44925, Peroxin 19 shRNA Plasmid (h): sc-44925-SH and Peroxin 19 shRNA (h) Lentiviral Particles: sc-44925-V.


Molecular Weight of Peroxin 19: 33 kDa.

Positive Controls: NTERA-2 cl.D1 whole cell lysate: sc-364181, U-698-M whole cell lysate: sc-364799 or T-47D cell lysate: sc-2293.

## RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

## DATA



Peroxin 19 (2E4): sc-517578. Western blot analysis of Peroxin 19 expression in NTERA-2 cl.D1 (A), U-698-M (B), T-47D (C) and PC-3 (D) whole cell lysates.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.