

HLA-F (201A2H): sc-517640

BACKGROUND

Major histocompatibility complex (MHC) molecules form an integral part of the immune response system. They are cell-surface receptors that bind peptides and present them to T lymphocytes. Human leukocyte antigens (HLAs) are polymorphic members of the MHC family that are specifically involved in the presentation of antigens to the T cell receptor. There are two classes of HLA antigens: class I and class II. Class I molecules are expressed in nearly all cells and play a central role in the immune system by presenting peptides derived from the endoplasmic reticulum. The differential structural properties of MHC class I and class II molecules account for their respective roles in activating different populations of T lymphocytes. HLA-F (major histocompatibility complex, class I, F), also known as CDA12, HLA-5.4 or HLA-CDA12, is a 346 amino acid single-pass type I membrane protein that belongs to the MHC class I family.

REFERENCES

1. Lury, D., Epstein, H. and Holmes, N. 1990. The human class I MHC gene HLA-F is expressed in lymphocytes. *Int. Immunol.* 2: 531-537.
2. Geraghty, D.E., Wei, X.H., Orr, H.T. and Koller, B.H. 1990. Human leukocyte antigen F (HLA-F). An expressed HLA gene composed of a class I coding sequence linked to a novel transcribed repetitive element. *J. Exp. Med.* 171: 1-18.
3. He, X., Xu, L., Liu, Y. and Zeng, Y. 2004. Identification of a novel HLA-F allele-HLA-F*010102. *Tissue Antigens* 63: 181-183.
4. Pyo, C.W., Williams, L.M., Moore, Y., Hyodo, H., Li, S.S., Zhao, L.P., Sageshima, N., Ishitani, A. and Geraghty, D.E. 2006. HLA-E, HLA-F, and HLA-G polymorphism: genomic sequence defines haplotype structure and variation spanning the nonclassical class I genes. *Immunogenetics* 58: 241-251.
5. Hasan, Z., Carlson, J.M., Gatanaga, H., Le, A.Q., Brumme, C.J., Oka, S., Brumme, Z.L. and Ueno, T. 2012. Minor contribution of HLA class I-associated selective pressure to the variability of HIV-1 accessory protein Vpu. *Biochem. Biophys. Res. Commun.* 421: 291-295.
6. Zeng, B., Wang, Y., Zhang, Y. and Liao, A.H. 2012. Mifepristone inhibits the progesterone-induced expressions of HLA-G, -E, -F genes in trophoblasts during first trimester. *Zhonghua Yi Xue Za Zhi* 92: 15-17.
7. Zhang, X., Lin, A., Zhang, J.G., Bao, W.G., Xu, D.P., Ruan, Y.Y. and Yan, W.H. 2013. Alteration of HLA-F and HLA I antigen expression in the tumor is associated with survival in patients with esophageal squamous cell carcinoma. *Int. J. Cancer* 132: 82-89.

CHROMOSOMAL LOCATION

Genetic locus: HLA-F (human) mapping to 6p22.1.

SOURCE

HLA-F (201A2H) is a mouse monoclonal antibody raised against recombinant HLA-F of human origin.

PRODUCT

Each vial contains 100 µg IgG₁ kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

HLA-F (201A2H) is recommended for detection of HLA-F of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000).

Suitable for use as control antibody for HLA-F siRNA (h): sc-105528, HLA-F shRNA Plasmid (h): sc-105528-SH and HLA-F shRNA (h) Lentiviral Particles: sc-105528-V.

Molecular Weight of HLA-F: 39 kDa.

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.