

Ricin A (RA519): sc-52189

BACKGROUND

Ricin, a type II ribosomal inactivating protein, inhibits protein biosynthesis by its RNA N-glycosidase activity. Ricin toxin, derived from the castor bean *Ricinus communis*, is a prototypic A-B toxin in which the B chain binds to the target cell, and the A chain (RTA) mediates the toxic activity. Ricin B chain (RTB) is a lectin that is responsible for cell agglutination and binds to β -D-galactopyranoside moieties found at the cell surface (e.g., on glycoproteins), allowing the A chain to enter the cell. In turn, the A chain functions enzymatically as an RNA N-glycosidase that depurinates adenine 4324 in the 28S rRNA of the 60S ribosomal subunit. The crystal structure of ricin has been defined.

REFERENCES

- Smith, D.C., Salio, M., Lord, J.M., Roberts, L.M. and Cerundolo, V. 2004. Lack of dendritic cell maturation by the plant toxin Ricin. *Eur. J. Immunol.* 34: 2149-2157.
- Spooner, R.A, Watson, P.D., Marsden, C.J., Smith, D.C., Moore, K.A., Cook, J.P., Lord, J.M. and Roberts, L.M. 2004. Protein disulphide-isomerase reduces Ricin to its A and B chains in the endoplasmic reticulum. *Biochem. J.* 383: 285-293.
- Moisenovich, M., Tonevitsky, A., Maljuchenko, N., Kozlovskaya, N., Agapov, I., Volkmandt, W. and Bereiter-Hahn, J. 2004. Endosomal Ricin transport: involvement of Rab4- and Rab5-positive compartments. *Histochem. Cell Biol.* 121: 429-439.
- Mantis, N.J., Farrant, S.A. and Mehta, S. 2004. Oligosaccharide side chains on human secretory IgA serve as receptors for Ricin. *J. Immunol.* 172: 6838-6845.
- Maddaloni, M., Cooke, C., Wilkinson, R., Stout, A.V., Eng, L. and Pincus, S.H. 2004. Immunological characteristics associated with the protective efficacy of antibodies to Ricin. *J. Immunol.* 172: 6221-6228.
- Amukele, T.K. and Schramm, V.L. 2004. Ricin A chain substrate specificity in RNA, DNA, and hybrid stem-loop structures. *Biochemistry* 43: 4913-4922.
- Bellisola, G., Fracasso, G., Ippoliti, R., Menestrina, G., Rosen, A., Solda, S., Udali, S., Tomazzolli, R., Tridente, G. and Colombatti, M. 2004. Reductive activation of Ricin and Ricin A chain immunotoxins by protein disulfide isomerase and thioredoxin reductase. *Biochem. Pharmacol.* 67: 1721-1731.
- Wu, Y.H., Shih, S.F. and Lin, J.Y. 2004. Ricin triggers apoptotic morphological changes through caspase-3 cleavage of BAT3. *J. Biol. Chem.* 279: 19264-19275.

SOURCE

Ricin A (RA519) is a mouse monoclonal antibody raised against agglutinin of *Ricinus communis* origin.

PRODUCT

Each vial contains 100 μ g IgG₁ in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

Ricin A (RA519) is recommended for detection of RCA120 and the A-chain of RCA60 of *Ricinus communis* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000).

Molecular Weight of Ricin A: 32 kDa.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.