

# AP-2 $\alpha$ / $\beta$ (SPM231): sc-52983

## BACKGROUND

AP-2 transcription factor family members include AP-2 $\alpha$ , AP-2 $\beta$  and AP-2 $\gamma$ , which specifically bind to the DNA consensus sequence CCCAGGC and initiate transcription of selected genes. AP-2, also known as ERF-1, plays a role in regulating estrogen receptor expression. AP-2 $\beta$ , a splice variant of AP-2 $\alpha$ , inhibits AP-2 activity. Besides subscribing to the AP-2 complex, AP-2 $\alpha$ , AP-2 $\beta$  and AP-2 $\gamma$  proteins compose the OB2-1 transcription factor complex. OB2-1 specifically upregulates expression of the proto-oncogene c-ErbB-2, which is overexpressed in 25-30% of breast cancers. The gene encoding AP-2 $\alpha$  maps to human chromosome 6p24. AP-2 $\alpha$  may play an important role in the development of ectodermal-derived tissues. Deleterious mutations involving the AP-2 $\alpha$  gene are linked to microphthalmia, corneal clouding and other anterior eye chamber defects. The ubiquitously expressed AP-4 transcription factor specifically binds to the DNA consensus sequence 5'-CAGCTG-3'. AP-4 interacts with promoters for immunoglobulin- $\kappa$  gene families and simian virus 40. AP-4 may enhance the transcription of the human Huntington's disease gene. AP-4 is a helix-loop-helix protein that contains two distinctive leucine repeat elements.

## REFERENCES

- Williams, T., Admon, A., Luscher, B. and Tjian, R. 1988. Cloning and expression of AP-2, a cell-type-specific transcription factor that activates inducible enhancer elements. *Genes Dev.* 2: 1557-1569.
- Buettner, R., Kannan, P., Imhof, A., Bauer, R., Yim, S.O., Glockshuber, R., Van Dyke, M.W. and Tainsky, M.A. 1993. An alternatively spliced mRNA from the AP-2 gene encodes a negative regulator of transcriptional activation by AP-2. *Mol. Cell. Biol.* 13: 4174-4185.
- Moser, M., Imhof, A., Pscherer, A., Bauer, R., Amselgruber, W., Sinowatz, F., Hofstadter, F., Schule, R. and Buettner, R. 1995. Cloning and characterization of a second AP-2 transcription factor: AP-2 $\beta$ . *Development* 121: 2779-2788.
- Bosher, J.M., Totty, N.F., Hsuan, J.J., Williams, T. and Hurst, H.C. 1996. A family of AP-2 proteins regulates c-ErbB-2 expression in mammary carcinoma. *Oncogene* 13: 1701-1707.
- Williamson, J.A., Bosher, J.M., Skinner, A., Sheer, D., Williams, T. and Hurst, H.C. 1996. Chromosomal mapping of the human and mouse homologues of two new members of the AP-2 family of transcription factors. *Genomics* 35: 262-264.
- McPherson, L.A., Baichwal, V.R. and Weigel, R.J. 1997. Identification of ERF-1 as a member of the AP2 transcription factor family. *Proc. Natl. Acad. Sci. USA* 94: 4342-4347.
- Moser, M., Ruschoff, J. and Buettner, R. 1997. Comparative analysis of AP-2 $\alpha$  and AP-2 $\beta$  gene expression during murine embryogenesis. *Dev. Dyn.* 208: 115-124.
- Satoda, M., Zhao, F., Diaz, G.A., Burn, J., Goodship, J., Davidson, H., Pierpont, M., Gelb, B. 2000. Mutations in TFAP2B cause Char syndrome, a familial form of patent ductus arteriosus. *Nature Genet.* 25: 42-46.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## CHROMOSOMAL LOCATION

Genetic locus: TFAP2A (human) mapping to 6p24.3.

## SOURCE

AP-2 $\alpha$ / $\beta$  (SPM231) is a mouse monoclonal antibody raised against a peptide from the C-terminus of AP-2 $\alpha$ .

## PRODUCT

Each vial contains 200  $\mu$ g IgG<sub>1</sub> in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

AP-2 $\alpha$ / $\beta$  (SPM231) is recommended for detection of AP-2 $\alpha$ / $\beta$  of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence and immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500).

Molecular Weight of AP-2 $\alpha$ : 48 kDa.

Molecular Weight of AP-2 $\beta$ : 47 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-mouse IgG-HRP: sc-2005 (dilution range: 1:2000-1:32,000) or Cruz Marker™ compatible goat anti-mouse IgG-HRP: sc-2031 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-mouse IgG-FITC: sc-2010 (dilution range: 1:100-1:400) or goat anti-mouse IgG-TR: sc-2781 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941. 3) Immunohistochemistry: use ImmunoCruz™: sc-2050 or ABC: sc-2017 mouse IgG Staining Systems.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.