



## TRF1 (TRF1 3H11): sc-53458

### BACKGROUND

Telomeric repeat binding factor 1 (TRF1, TERF1, PIN2, TRBF1) and telomeric repeat binding factor 2 (TRF2, TERF2, TRBF2) are present at telomeres throughout the cell cycle, where they regulate telomerase by acting in *cis* to limit the elongation of individual chromosome ends. Telomerase adds hexameric repeats of 5'-TTAGGG-3' to the ends of chromosomal DNA. This telomerase enzyme plays an influential role in cellular immortalization and cellular senescence. TRF1 negatively regulates telomere elongation, while TRF2 protects the chromosome ends by inhibiting end-to-end fusions. Downregulation of TRF expression in tumor cells may contribute to cell immortalization and malignant progression. TRF1 has an acidic N-terminus while TRF2 has a basic N-terminus. TRF2 localizes in the nucleolus at G<sub>0</sub> and S and diffuses out of the nucleolus in G<sub>2</sub> phase. During mitosis TRF2 disperses from the condensed chromosomes and returns to the nucleolus at cytokinesis.

### REFERENCES

1. Aragona, M., Buda, C.A., Panetta, S., Morelli, M., Giudice, A., Campagna, F.L., Pontoriero, A., Cascinu, S. and La Torre, F. 2000. Immunohistochemical telomeric-repeat binding factor-1 expression in gastrointestinal tumors. *Oncol. Rep.* 7: 987-990.
2. Matsutani, N., Yokozaki, H., Tahara, E., Tahara, H., Kuniyasu, H., Haruma, K., Chayama, K., Yasui, W. and Tahara, E. 2001. Expression of telomeric repeat binding factor 1 and 2 and TRF1-interacting nuclear protein 2 in human gastric carcinomas. *Int. J. Oncol.* 19: 507-512.
3. Yajima, T., Yagihashi, A., Kameshima, H., Kobayashi, D., Hirata, K. and Watanabe, N. 2001. Telomerase reverse transcriptase and telomeric-repeat binding factor protein 1 as regulators of telomerase activity in pancreatic cancer cells. *Br. J. Cancer* 85: 752-757.
4. Seimiya, H. and Smith, S. 2002. The telomeric poly(ADP-ribose) polymerase, tankyrase 1, contains multiple binding sites for telomeric repeat binding factor 1 (TRF1) and a novel acceptor, 182 kDa tankyrase-binding protein (TAB182). *J. Biol. Chem.* 277: 14116-14126.
5. Nakanishi, K., Kawai, T., Kumaki, F., Hiroi, S., Mukai, M., Ikeda, E., Koering, C.E. and Gilson, E. 2003. Expression of mRNAs for telomeric repeat binding factor (TRF)-1 and TRF2 in atypical adenomatous hyperplasia and adenocarcinoma of the lung. *Clin. Cancer Res.* 9: 1105-1111.
6. Yang, S.W., Kim, D.H., Lee, J.J., Chun, Y.J., Lee, J.H., Kim, Y.J., Chung, I.K. and Kim, W.T. 2003. Expression of the telomeric repeat binding factor gene NgTRF1 is closely coordinated with the cell division program in tobacco BY-2 suspension culture cells. *J. Biol. Chem.* 278: 21395-21407.
7. Zhang, S., Hemmerich, P. and Grosse, F. 2004. Nucleolar localization of the human telomeric repeat binding factor 2 (TRF2). *J. Cell. Sci.* 117: 3935-3945.
8. La Torre, D., de Divitiis, O., Conti, A., Angileri, F.F., Cardali, S., Aguenouz, M., Aragona, M., Panetta, S., d'Avella D, Vita G, La Torre F, Tomasello F. 2005. Expression of telomeric repeat binding factor-1 in astroglial brain tumors. *Neurosurgery* 56: 802-810.

### CHROMOSOMAL LOCATION

Genetic locus: TERF1 (human) mapping to 8q13; Terf1 (mouse) mapping to 1 E3.

### SOURCE

TRF1 (TRF1 3H11) is a mouse monoclonal antibody raised against human TRF1.

### PRODUCT

Each vial contains 200 µg IgG<sub>1</sub> in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

### APPLICATIONS

TRF1 (TRF1 3H11) is recommended for detection of TRF1 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

Suitable for use as control antibody for TRF1 siRNA (h): sc-36722.

Molecular Weight of TRF1: 60 kDa.

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-mouse IgG-HRP: sc-2005 (dilution range: 1:2000-1:32,000) or Cruz Marker™ compatible goat anti-mouse IgG-HRP: sc-2031 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-mouse IgG-FITC: sc-2010 (dilution range: 1:100-1:400) or goat anti-mouse IgG-TR: sc-2781 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

### STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.