

IL-10R (1B1.3a): sc-53706

BACKGROUND

The IL-10 receptor, IL-10R, is a member of the class II subgroup of the cytokine receptor family and exhibits structural similarity to the interferon receptor. IL-10R is expressed in B cells and T helper cells, as well as in LPS-induced mouse fibroblasts. Overall, mouse IL-10R and human IL-10R share 60% sequence identity at the protein level. Stimulation with IL-10 leads to phosphorylation of JAK1 and Tyk 2 tyrosine kinases. The activated kinases phosphorylate the two tyrosine residues (Tyr 446 and Tyr 496) in the cytoplasmic domain of IL-10R α . The phosphorylation of these two residues are required for proper function of IL-10R and activation of IL-10E1 signaling. IL-10R β is ubiquitously expressed and, in addition to forming the IL-10 heterodimeric receptor, it forms a heterodimeric receptor with an IL-22R subunit and an IL-28R subunit. IL-10R is constitutively expressed on human natural killer (NK) cells and the direct binding of IL-10 potentiates cytokine production by human NK cells.

REFERENCES

1. Rousset, F., Garcia, E., Defrance, T., Peronne, C., Vezzio, N., Hsu, D.H., Kastelein, R., Moore, K.W. and Banchereau, J. 1992. Interleukin-10 is a potent growth and differentiation factor for activated human B lymphocytes. *Proc. Natl. Acad. Sci. USA* 89: 1890-1893.
2. Ho, A.S.Y., Liu, Y., Khan, T.A., Hsu, D.H., Bazan, J.F. and Moore, K.W. 1993. A receptor for interleukin-10 is related to interferon receptors. *Proc. Natl. Acad. Sci. USA* 90: 11267-11271.
3. Liu, Y., Wei, S.H., Ho, A.S., de Waal Malefyt, R. and Moore, K.W. 1994. Expression cloning and characterization of a human IL-10 receptor. *J. Immunol.* 152: 1821-1829.
4. Weber-Nordt, R.M., Meraz, M.A. and Schreiber, R.D. 1994. Lipopolysaccharide-dependent induction of IL-10 receptor expression on murine fibroblasts. *J. Immunol.* 153: 3734-3744.
5. Ho, A.S., Wei, S.H., Mui, A.L., Miyajima, A. and Moore, K.W. 1995. Functional regions of the mouse interleukin-10 receptor cytoplasmic domain. *Mol. Cell. Biol.* 15: 5043-5053.
6. Tan, J.C., Braun, S., Rong, H., DiGiacomo, R., Dolphin, E., Baldwin, S., Narula, S.K., Zavodny, P.J. and Chou, C.C. 1995. Characterization of recombinant extracellular domain of human interleukin-10 receptor. *J. Biol. Chem.* 270: 12906-12911.
7. Carson, W.E., Lindemann, M.J., Baiocchi, R., Linett, M., Tan, J.C., Chou, C.C., Narula, S. and Caligiuri, M.A. 1995. The functional characterization of interleukin-10 receptor expression on human natural killer cells. *Blood* 85: 3577-3585.

CHROMOSOMAL LOCATION

Genetic locus: Il10rb (mouse) mapping to 16 C3.3.

SOURCE

IL-10R (1B1.3a) is a rat monoclonal antibody raised against soluble IL-10R of mouse origin.

PRODUCT


Each vial contains 100 μ g IgG₁ in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

IL-10R (1B1.3a) is recommended for detection of IL-10R of mouse origin by flow cytometry (1 μ g per 1×10^6 cells).

Molecular Weight of IL-10R: 90-110 kDa.

DATA



IL-10R (1B1.3a): sc-53706. Indirect FCM analysis of C57BL/6 splenocytes stained with IL-10R (1B1.3a), followed by PE-conjugated goat anti-rat IgG: sc-3740. Black line histogram represents the isotype control, normal rat IgG₁: sc-3882.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.