



KLH (MTK2113): sc-53715

BACKGROUND

Hemocyanin proteins, which are freely dissolved in mollusk and arthropod hemolymph, belong to the tyrosinase family of proteins. Keyhole limpet hemocyanin (KLH) is an extracellular respiratory protein that is commonly used as hapten carrier. It is also used as an immune stimulant since it activates T cells and induces higher titres of IgE and IgG₁ antibody levels. The marine gastropod *Megathura crenulata* expresses KLH. KLH has two isoforms: KLH1 and KLH2 of 400 kDa and 345 kDa, respectively. KLH is useful as an immunotherapeutic alternative for patients with cancer, and it also has implications in the treatment of drug addiction.

REFERENCES

- Harris, J.R. and Markl, J. 1999. Keyhole limpet hemocyanin (KLH): a biomedical review. *Micron* 30: 597-623.
- Holmberg, L.A., Oparin, D.V., Gooley, T., Lilleby, K., Bensinger, W., Reddish, M.A., MacLean, G.D., Longenecker, B.M, and Sandmaier, B.M. 2000. Clinical outcome of bre chemotherapy, autologous stem cell rescue and THERATOPE STn-KLH cancer vaccine. *Bone Marrow Transplant*. 25: 1233-1241.
- Schumacher, K. 2001. Keyhole limpet hemocyanin (KLH) conjugate vaccines as novel therapeutic tools in malignant disorders. *J. Cancer Res. Clin. Oncol.* 127: R1-R2.
- Musselli, C., Livingston, P.O. and Ragupathi, G. 2001. Keyhole limpet hemocyanin conjugate vaccines against cancer: the Memorial Sloan Kettering experience. *J. Cancer Res. Clin. Oncol.* 127 : R20-R206.
- Gebauer, W., Robin Harris, J. and Markl, J. 2002. Topology of the 10 subunits within the decamer of KLH, the hemocyanin of the marine gastropod *Megathura crenulata*. *J. Struct. Biol.* 139: 153-159.
- Höglund, A.S., Jones, A.M. and Josefsson, L.G. 2002. An antigen expressed during plant vascular development crossreacts with antibodies towards KLH (keyhole limpet hemocyanin). *J. Histochem. Cytochem.* 50: 999-1003.
- Kurokawa, T., Wuhrer, M., Lochnit, G., Geyer, H., Markl, J. and Geyer, R. 2002. Hemocyanin from the keyhole limpet *Megathura crenulata* (KLH) carries a novel type of N-glycans with Gal(beta1-6) Man-motifs. *Eur. J. Biochem.* 269: 5459-5473.
- Millard, A.L., Ittelet, D., Schooneman, F. and Bernard, J. 2003. Dendritic cell KLH loading requirements for efficient CD4⁺ T cell priming and help to peptide-specific cytotoxic T cell response, in view of potential use in cancer vaccines. *Vaccine* 21: 869-876.
- Roman, D., Ulrich, P., Paul, G., Court, M., Vit, P., Kehren, J. and Mahl, A. 2004. Determination of the effect of calcineurin inhibitors on the rat's immune system after KLH immunisation. *Toxicol. Lett.* 149: 133-140.

SOURCE

KLH (MTK2113) is a mouse monoclonal antibody raised against KLH.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 100 µg IgG₁ in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

KLH (MTK2113) is recommended for detection of KLH by solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.