

LAT1 (Q-20): sc-54229

BACKGROUND

L-type amino acid transporter 1 (LAT1) is a multipass-membrane protein responsible for sodium-independent, high-affinity transport of large neutral amino acids. LAT1 functions as a disulfide-linked heterodimer with the amino acid transport protein CD98. LAT1 is expressed predominantly in adult lung and liver but is also expressed in brain, thymus, retina, testis, placenta, bone marrow and fetal liver. In the retina, LAT1 localizes to the blood-retinal-barrier (BRB) and mediates L-leucine transport from the blood to the retina. The devastating effects on the brain caused by phenylketonuria are due to the increased levels of LAT1 on the blood-brain-barrier in response to high concentrations of phenylalanine in the blood. LAT1 accepts the amino-acid related anticancer agent melphalan and plays a significant role in cell proliferation, differentiation, and invasion in esophageal squamous cell carcinoma.

REFERENCES

1. Satoh, S., et al. 1997. Involvement of L-type-like amino acid transporters in S-nitrosocysteine-stimulated noradrenaline release in the rat hippocampus. *J. Neurochem.* 69: 2197-2205.
2. Yanagida, O., et al. 2001. Human L-type amino acid transporter 1 (LAT1): characterization of function and expression in tumor cell lines. *Biophys. Acta* 1514: 291-302.
3. Tamai, S., et al. 2001. Expression of L-type amino acid transporter 1 in a rat model of liver metastasis: positive correlation with tumor size. *Cancer Detect. Prev.* 25: 439-445.
4. Ohkame, H., et al. 2001. Expression of L-type amino acid transporter 1 (LAT1) and 4F2 heavy chain (4F2hc) in liver tumor lesions of rat models. *J. Surg. Oncol.* 78: 265-271.
5. Uchino, H., et al. 2002. Transport of amino acid-related compounds mediated by L-type amino acid transporter 1 (LAT1): insights into the mechanisms of substrate recognition. *Mol. Pharmacol.* 61: 729-737.
6. Boado, R.J., et al. 2004. Human LAT1 single nucleotide polymorphism N230K does not alter phenylalanine transport. *Mol. Genet. Metab.* 83: 306-311.
7. Tomi, M., et al. 2005. L-type amino acid transporter 1-mediated L-leucine transport at the inner blood-retinal barrier. *Invest. Ophthalmol. Vis. Sci.* 46: 2522-2530.
8. Kobayashi, H., et al. 2005. Expression of L-type amino acid transporter 1 (LAT1) in esophageal carcinoma. *J. Surg. Oncol.* 90: 233-238.
9. Li, S. and Whorton, A.R. 2005. Identification of stereoselective transporters for S-nitroso-L-cysteine: role of LAT1 and LAT2 in biological activity of S-nitrosothiols. *J. Biol. Chem.* 280: 20102-20110.

CHROMOSOMAL LOCATION

Genetic locus: Slc7a5 (mouse) mapping to 8 E1.

SOURCE

LAT1 (Q-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of LAT1 of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-54229 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

LAT1 (Q-20) is recommended for detection of L-type amino acid transporter 1 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for LAT1 siRNA (m): sc-62556, LAT1 shRNA Plasmid (m): sc-62556-SH and LAT1 shRNA (m) Lentiviral Particles: sc-62556-V.

Molecular Weight (observed) of LAT1: 40 kDa.

Molecular Weight of glycosylated LAT1: 45-55 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.