

Histone H2A.Z/H2A.F/Z (G-13): sc-54387

BACKGROUND

Histone H2A.Z/H2A.F/Z (H2A/z) is a 128 amino acid protein encoded by the human gene H2AFZ. Eukaryotic histones are basic and water soluble nuclear proteins that form hetero-octameric nucleosome particles by wrapping 146 base pairs of DNA sequentially in a left-handed super-helical turn to form chromosomal fiber. Two molecules of each of the four core histones (H2A, H2B, H3 and H4) form the octamer, which is comprised of two H2A-H2B dimers and two H3-H4 dimers, creating two nearly symmetrical halves by tertiary structure. H2A.Z/H2A.F/Z is a variant Histone H2A which replaces conventional H2A in a subset of nucleosomes. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of posttranslational modifications of histones, also called histone code, and nucleosome remodeling. H2A.Z/H2A.F/Z may be involved in the formation of constitutive heterochromatin and may be required for chromosome segregation during cell division.

REFERENCES

- Martin, C., et al. 2005. The diverse functions of histone lysine methylation. *Nat. Rev. Mol. Cell Biol.* 6: 838-849.
- Bustin, M., et al. 2005. The dynamics of Histone H1 function in chromatin. *Mol. Cell* 17: 617-620.
- de la Cruz, X., et al. 2005. Do protein motifs read the histone code? *Bioessays* 27: 164-175.
- Gunjan, A., et al. 2005. Regulation of histone synthesis and nucleosome assembly. *Biochimie* 87: 625-635.
- Rupp, R.A., et al. 2005. Gene regulation by Histone H1: new links to DNA methylation. *Cell* 123: 1178-1179.
- Bode, A.M., et al. 2005. Inducible covalent posttranslational modification of Histone H3. *Sci. STKE* 2005: re4.
- Nightingale, K.P., et al. 2006. Histone modifications: signaling receptors and potential elements of a heritable epigenetic code. *Curr. Opin. Genet. Dev.* 16: 125-136.

CHROMOSOMAL LOCATION

Genetic locus: H2AFZ (human) mapping to 4q23, H2afz (human) mapping to 7p13; H2afz (mouse) mapping to 3 G3, H2AFZ (mouse) mapping to 11 A1.

SOURCE

Histone H2A.Z/H2A.F/Z (G-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of histone H2A.Z of mouse origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-54387 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-54387 X, 200 µg/0.1 ml.

APPLICATIONS

Histone H2A.Z/H2A.F/Z (G-13) is recommended for detection of histone H2A.Z and histone H2A.F/Z of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Histone H2A.Z/H2A.F/Z (G-13) is also recommended for detection of histone H2A.Z and histone H2A.F/Z in additional species, including equine, canine, bovine, porcine and avian.

Histone H2A.Z/H2A.F/Z (G-13) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of Histone H2A.Z/H2A.F/Z: 14 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200, mouse placenta extract: sc-364247 or rat placenta extract: sc-364808.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.