

# CA VB (T-17): sc-54764

## BACKGROUND

Carbonic anhydrases (CAs) are members of a large family of zinc metalloenzymes responsible for catalyzing the reversible hydration of carbon dioxide. CAs show extensive diversity in their distribution and subcellular localization. They are involved in a variety of biological processes, including calcification, bone resorption, respiration, acid-base balance and the formation of aqueous humor, saliva, gastric juice and cerebrospinal fluid. CA VB, also known as carbonate dehydratase VB, is one of two isoforms of CA V. It localizes to the mitochondria and is involved in metabolic processes. CA VB is predominantly expressed in heart, pancreas, lung, placenta, kidney and skeletal muscle. It exhibits highest homology with family member CA VA (the second isoform of CA V); however, unlike CA VA, it is not expressed in the liver, suggesting that it plays a significantly different physiological role.

## REFERENCES

1. Fujikawa-Adachi, K., et al. 1999. Human mitochondrial carbonic anhydrase VB. cDNA cloning, mRNA expression, subcellular localization and mapping to chromosome X. *J. Biol. Chem.* 274: 21228-21233.
2. Shah, G.N., et al. 2000. Mitochondrial carbonic anhydrase CA VB: differences in tissue distribution and pattern of evolution from those of CA VA suggest distinct physiological roles. *Proc. Natl. Acad. Sci. USA* 97: 1677-1682.
3. Nishimori, I., et al. 2001. Carbonic anhydrase isozymes in the human pancreas. *Dig. Liver Dis.* 33: 68-74.
4. Winum, J.Y., et al. 2006. Carbonic anhydrase inhibitors: clash with Ala 65 as a means for designing inhibitors with low affinity for the ubiquitous isozyme II, exemplified by the crystal structure of the topiramate sulfamide analogue. *J. Med. Chem.* 49: 7024-7031.
5. Supuran, C.T. 2007. Carbonic anhydrases as drug targets—an overview. *Curr. Top. Med. Chem.* 7: 825-833.
6. De Simone, G., et al. 2007. Antiobesity carbonic anhydrase inhibitors. *Curr. Top. Med. Chem.* 7: 879-884.
7. Vitale, R.M., et al. 2007. Molecular modeling study for the binding of zonisamide and topiramate to the human mitochondrial carbonic anhydrase isoform VA. *Bioorg. Med. Chem.* 15: 4152-4158.
8. Supuran, C.T., et al. 2007. Carbonic anhydrases as targets for medicinal chemistry. *Bioorg. Med. Chem.* 15: 4336-4350.
9. Temperini, C., et al. 2007. Carbonic anhydrase inhibitors. Interaction of the antiepileptic drug sultiame with twelve mammalian isoforms: Kinetic and X-ray crystallographic studies. *Bioorg. Med. Chem. Lett.* 17: 4866-4872.

## CHROMOSOMAL LOCATION

Genetic locus: CA5B (human) mapping to Xp22.2.

## SOURCE

CA VB (T-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of CA VB of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-54764 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

CA VB (T-17) is recommended for detection of CA VB of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

CA VB (T-17) is also recommended for detection of CA VB in additional species, including equine, canine and porcine.


Suitable for use as control antibody for CA VB siRNA (h): sc-62034, CA VB shRNA Plasmid (h): sc-62034-SH and CA VB shRNA (h) Lentiviral Particles: sc-62034-V.

Molecular Weight of CA VB precursor: 36 kDa.

Molecular Weight of mature CA VB: 32 kDa.

Positive Controls: A-673 cell lysate: sc-2414 or HEK293 whole cell lysate: sc-45136.

## DATA



CA VB (T-17): sc-54764. Western blot analysis of CA VB expression in A-673 (A) and HEK293 (B) whole cell lysates.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.